**Cordis**<sup>TM</sup>

100000002391.5

Instructions for Use

**SABER<sup>TM</sup> Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter**

Mode d'emploi

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
## English

**STERILE. Sterilized with ethylene oxide gas. Nonpyrogenic. Radiopaque. For single use only. Do not resterilize.**

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

### I. Device name

The device brand name is **SABER™** Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter; the generic device name is PTA Balloon Dilatation Catheter.

### II. Description

The Cordis **SABER™** Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter is a catheter with a distal inflatable balloon. Radiopaque marker bands indicate the dilating section of the balloon and aid in balloon placement. For balloon lengths greater than or equal to 100mm, the distal section will have one (1) marker band and proximal section will consist of two (2) adjacent marker bands. For balloon lengths less than 100mm, the distal and proximal section will have each one (1) marker band. The catheter tip is tapered to ease entry into peripheral arteries and to facilitate the crossing of tight stenoses. **SABER™** balloons are coated with a hydrophilic material designed to increase lubricity throughout the lifetime of the device.

The working pressure range for the balloon is between the nominal pressure and the rated burst pressure. All balloons distend to sizes above the nominal size at pressures greater than the nominal pressure. Consult the compliance table provided below, and also incorporated with the product for typical balloon diameters at given pressures.

**Compliance Table. Balloon Diameter vs. Inflation Pressure**

Inflation Pressure		Ø Balloon Diameter (mm)										
atm	KPa	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm	9.0 mm	10.0 mm
6	608	1.87	2.35	2.91	3.35	3.74	4.67	5.75	6.86	7.73	8.81	9.87
7	709	1.91	2.40	2.97	3.43	3.83	4.79	5.89	6.99	7.88	8.96	10.02
8	811	1.95	2.44	3.01	3.49	3.91	4.88	5.98	7.08	8.00	9.07	10.14
9	912	1.97	2.47	3.04	3.54	3.98	4.96	6.06	7.16	8.10	9.17	10.24
10	1013	2.00	2.51	3.08	3.58	4.04	5.02	6.13	7.24	8.18	9.25	10.34
11	1115	2.03	2.53	3.10	3.62	4.08	5.07	6.19	7.31	8.26		
12	1216	2.06	2.56	3.13	3.66	4.13	5.13	6.25	7.37	8.33		
13	1317	2.09	2.58	3.16	3.68	4.17	5.17	6.31	7.43			
14	1419	2.11	2.61	3.18	3.72	4.20	5.22	6.37	7.50			
15	1520	2.13	2.63	3.21	3.75	4.24	5.26	6.43				
16	1621	2.15	2.65	3.23	3.78	4.27	5.31	6.49				
17	1723	2.18	2.67	3.25	3.81	4.31						
18	1824	2.20	2.70	3.28	3.84	4.34						

Nom.P  
Nominal Pressure

RBP  
RBP  
Rated Burst Pressure

RBP  
RBP  
RBP Rated Burst Pressure  
(For balloon length of 150-300mm)

**NOTE: Do not exceed the RBP.**

The balloon lumen, marked on the hub as "BALLOON" is used to inflate and deflate the balloon. The nominal balloon size is printed on the hub. The guidewire lumen, marked "THRU" is used to track the catheter over a prepositioned guidewire.

The radiopaque marker bands (the length between the outer edge of most proximal to the outer edge of most distal marker band) indicate the stated nominal length of the balloon.

### III. Indications

The **SABER™** PTA catheter is intended to dilate stenoses in iliac, femoral, ilio-femoral, popliteal and infra-popliteal arteries. The device is also indicated for post-dilation of balloon-expandable and self-expanding stents in the peripheral vasculature.

### IV. Contraindications

- None known for PTA procedure.
- The **SABER™** PTA Balloon Dilatation catheter is contraindicated for use in coronary arteries.

### V. Warnings

- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- Exposure to temperatures above 54°C (130°F) may damage the device.
- Do not expose the device to organic solvents (e.g. alcohol).
- Do not use Ethiodol or Lipiodol® contrast media.
- Store in a cool, dark, dry place.
- Use the catheter prior to the "Use By" date specified on the package.
- Do not use if inner package is opened or damaged.
- Prior to angioplasty, the catheter should be examined to verify functionality and integrity, and ensure that its size and shape are suitable for the specific procedure for which it is to be used. Do not use if product damage is suspected or evident.
- To reduce the potential for vessel damage or the risk of dislodgement of particles it is very important that the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the lesion. The balloon dimensions are printed on the product label. The compliance table incorporated with the product shows how balloon diameter increases as pressure increases.
- Do not exceed the rated burst pressure recommended on the label. The rated burst pressure is based on the results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure. Use of a pressure monitoring device is recommended to prevent over-pressurization.
- Pressure in excess of the rated burst pressure can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath. Balloon rupture can cause vessel damage and the need for additional intervention.
- Use only the recommended balloon inflation medium (a 50/50 mixture by volume of contrast medium and normal saline). Never use air or any gaseous medium to inflate the balloon.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of resistance before proceeding.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a

- unit, particularly if balloon rupture or leakage is known or suspected.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter. Forceful handling can result in catheter separation and the subsequent need to use a snare or other medical interventional techniques to retrieve the pieces.
- Always verify integrity of the catheter after removal.
- Maintain a snug seal with the hemostasis valve over the balloon catheter during advancement to prevent introduction of air into the sheath or guiding catheter. Without a snug seal, a tight fit between the balloon section of the balloon catheter and the sheath or guiding catheter may cause a risk for introduction of air and air entrapment during advancement of the balloon catheter through the sheath or guiding catheter.
- The minimal acceptable sheath introducer/guiding catheter size is printed on the package label. Do not attempt to pass the PTA catheter through smaller size sheath introducer/guiding catheter than indicated on the label. Use of a smaller than indicated accessory device can lead to introduction of air into that device as the balloon catheter is advanced, which may not be removed during air aspiration.
- Do not use with power injectors.

### VI. Precautions

- The catheter system should be used only by physicians trained in the performance of arteriography and who have received appropriate training in percutaneous transluminal angioplasty.
- Flush all devices entering the vascular system with sterile heparinized saline or similar isotonic solution. Prior to use, ensure all devices have been flushed and air is removed from the system according to standard medical practice. Failure to do so could result in air entering the vascular system.
- Caution should be taken when treating patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.
- During the procedure remove blood or any other residues from the devices, using heparinized-saline soaked gauze.
- The catheter system is not intended for precise arterial blood pressure monitoring.

### VII. Complications

Potential complications include but are not limited to:

- Abrupt vessel closure
- Additional intervention
- Acute myocardial infarction
- Allergic reaction (device, contrast medium and medications)
- Amputation
- Arrhythmias
- Arteriovenous fistula
- Bradycardia
- Death
- Embolism
- Hematoma at puncture site
- Hemorrhage
- Hypotension / hypertension
- Inflammation/ infection/ sepsis
- Ischemia
- Necrosis
- Neurological events , including peripheral nerve injury, transient ischemic attack and/ or stroke
- Organ failure (single, multiple)
- Pain
- Paralysis
- Potential for balloon burst and potential complications (rated burst pressure)
- Potential for separation and potential complications (integrity to be checked before and after use)
- Procedural complications: bleeding, hypotension, access site complications
- Pseudoaneurysm
- Renal failure
- Restenosis of the dilated vessel
- Thrombosis
- Vascular complications (e.g. intimal tear, dissection, pseudoaneurysm, perforation, rupture, spasm, occlusion)

### VIII. Instructions for Use

#### Removal from Package

Open the pouch, grasp the hub and gently take the catheter out.

#### Preparation

- Attach a 3-way stopcock to the inflation port, which is marked "BALLOON".
  - Attach a partially filled syringe with heparinized saline to the stopcock, open the stopcock to the balloon and induce negative pressure.
  - Hold the syringe and proximal end of the catheter above the distal end of the catheter, and hold the balloon vertically with the balloon tip pointing down.
  - While maintaining negative pressure close the stopcock to the inflation port. Remove the syringe and purge the air.
  - To ensure air contained in the balloon and inflation lumen is removed, apply negative pressure twice as instructed and repeat steps 2-4.
  - Without twisting, slide the forming tube off the balloon.
  - Prepare an angioplasty inflation system with a 50% solution of contrast medium in sterile saline or similar solution.
  - Purge the air from the inflation device.
  - Connect the inflation device to the 3-way stopcock that is connected to the catheter inflation port, open the stopcock to the catheter and slowly fill the inflation lumen and the balloon will slowly fill with diluted contrast medium.
- Caution:** Do not apply positive or negative pressure to the balloon at this time.

#### Insertion, Inflation and Withdrawal

- Flush the "THRU" lumen with sterile heparinized saline or a similar isotonic solution.
- Place the prepared catheter over a prepositioned guidewire and advance the tip to the introduction site.
  - Note:** Balloon inflation should be performed with the guidewire extended beyond the catheter tip. It is strongly recommended that the guidewire, the balloon catheter, or both, remain across the lesion until the procedure is complete and the dilatation system is to be removed from the vessel.
  - Note:** To preserve the folded balloon shape during insertion and catheter manipulation, maintain a vacuum on the inflation lumen.
  - Caution:** Fully deflate the balloon by inducing negative pressure with the inflation system whenever the PTA catheter is advanced or withdrawn. Do not advance or withdraw the PTA catheter within the vasculature unless the catheter is preceded by a guidewire.
- Carefully advance the catheter through a sheath or guide catheter through the percutaneous entry site.
  - Note:** Gentle counterclockwise rotation of the balloon may ease introduction through the sheath or percutaneous entry site.
  - Note:** Perform all further catheter manipulations under fluoroscopy.
- Carefully advance the catheter to the selected stenosis.
  - Caution:** If strong resistance is met during advancement or withdrawal of the catheter, discontinue movement and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system.

\* Ethiodol and Lipiodol are trademarks of Guerbet S.A.

- Using fluoroscopy and the radiopaque marker bands, position the catheter at the appropriate location.
- When an acceptable position has been obtained, inflate the balloon to achieve the desired dilatation.  
**Caution:** Do not exceed the rated burst pressure. Higher pressures may damage the balloon or catheter or overdistend the selected artery.  
**Warning:** Inflation at a high rate may damage the balloon.
- Deflate the balloon by pulling vacuum on the inflation syringe or inflation device.
- Remove the vacuum (do not apply pressure) and carefully withdraw and remove the catheter.  
**Note:** Gentle counterclockwise rotation of the balloon may ease withdrawal from the sheath or from the percutaneous entry site. If the balloon cannot be withdrawn through the sheath, withdraw the catheter and sheath as a unit.

Dispose of all used devices in accordance with hospital policy for bio hazardous materials.

#### IX. References

The physician should consult recent literature on current medical practice on balloon dilatation.

#### X. Disclaimer of Warranty and Limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

THIS PAGE HAS BEEN INTENTIONALLY  
LEFT BLANK

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14201 NW 60th Avenue, Miami Lakes,  
Florida – Estados Unidos da América

Cordis Cashel  
Cashir Road, Cashel, County,  
Tipperary, Irlanda

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**CORDIS MEDICAL BRASIL LTDA.**

Rua Ponta Porã, nº 371, Vila Ipojuca

São Paulo - SP – 05.058-000

CNPJ: 27.548.227/0001-22

Farm. Resp: Alexandra Peixoto Demori  
Lima – CRF/SP: 42.520

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Patents: [www.cordis.com/na/patents](http://www.cordis.com/na/patents)

December 2022

10000002391.5

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Cordis US Corp. (também conhecida como