


# MYNXGRIP™

## VASCULAR CLOSURE DEVICE



### Ordering Information

The **MYNXGRIP™ VCD** includes:

- (1) Balloon catheter with integrated sealant
- (1) 10ml locking syringe

Size	Color	Order Number
6F/7F	Green	MX6721
5F	Gray	MX5021

1. FINALE Study Report: Clinical Results, TPR 1814-02. ‡ Time to discharge eligibility as compared to manual compression. MATRIX Clinical Trial (IDE# G030182).  
2. Hutchings D, Hayat A, et al. Success, Safety, and Efficacy of the Mynx Femoral Closure Device in a Real-World Cohort: Single-Center Experience. J Invasive Cardiol 2016 Mar; 28(3):104-8.  
3. Diamantopoulos A, Nourzaie R, et. al. Safety and efficacy of the Mynx Control vascular closure device in peripheral arterial procedures: A prospective study. Vascular. 2021; Dec 27; 17085381211062745. doi: 10.1177/17085381211062745

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO, MYNXGRIP and GRIP TECHNOLOGY are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2024 Cordis. All Rights Reserved. 100526838-3 07/2024



# MOVE FORWARD WITH MYNX




**MYNXGRIP™**  
VASCULAR CLOSURE DEVICE




# WHY COMPROMISE?

The MYNXGRIP™ Vascular Closure Device provides secure vascular closure without compromise.

The MYNXGRIP™ VCD offers the security of mechanical closure combined with the **safety of an extravascular sealant**. The MYNXGRIP™ VCD offers a patient-friendly closure option with no sutures, clamping, or metal implants and dissolves within 30 days leaving nothing behind but a healed artery or vein.



## Deployment Steps




Achieve temporary hemostasis and position at the arteriotomy or venotomy

The Grip Tip securely adheres to the artery or vein and the sealant fills the tissue tract


Platelets and blood cells collect inside the sealant's porous matrix

The sealant dissolves within 30 days leaving nothing behind but a healed artery or vein

# SEALANT SCIENCE FEATURING GRIP TECHNOLOGY




The sealant in the MYNXGRIP™ VCD consists of Polyethylene Glycol (PEG), a water-soluble, bio-inert, non-thrombogenic polymer, and is comprised of two configurations of PEG, the Grip Tip and the sealant.



Once the sealant enters the tissue tract, the body's temperature and pH level cause the Grip Tip to soften and securely adhere to the vessel wall, effectively gripping the artery and providing active closure.



The sealant's porous structure absorbs blood and subcutaneous fluids. The sealant swells three to four times its original size, filling the tissue tract.



The sealant actively adheres to the artery or vein while expanding and filling the tissue tract, providing a durable hemostasis and a platform for natural healing.

## SECURITY OF MECHANICAL CLOSURE, SAFETY OF AN EXTRAVASCULAR SEALANT.