


# DEPLOYMENT GUIDE




**MYNX CONTROL™ VENOUS**  
VASCULAR CLOSURE DEVICE



# SECURE CLOSURE RESULTS


POSITION THE BALLOON

**RESULT: TEMPORARY HEMOSTASIS**




DEPLOY THE SEALANT

**RESULT: CONSISTENT DELIVERY**



REMOVE THE DEVICE

**RESULT: SECURE EXTRAVASCULAR CLOSURE**



## PREPARATION

1



### PREPARE BALLOON

- Fill locking syringe with 2 to 3 mL of sterile saline, attach to stopcock and draw vacuum. Check luer connector and tighten if necessary.
- Inflate the balloon until the black marker on the inflation indicator is fully visible. Check for leaks in the balloon and syringe connector; retighten if necessary. Do not use the device if the balloon does not maintain pressure.
- Check for air bubbles in the balloon. If air bubbles are visible, deflate the balloon, draw vacuum to remove bubbles and re-inflate.
- Deflate balloon and leave syringe at neutral.

## POSITION THE BALLOON

2A



### INSERT DEVICE

- Insert the MYNX CONTROL™ VENOUS Vascular Closure Device into the procedural sheath through the sheath valve.
- Advance the catheter until the sheath catch nears the hub of the sheath.
- Rotate the sheath catch as needed to hook onto the side port of the procedural sheath.


2B



### INFLATE BALLOON

- Inflate the balloon until the black marker is fully visible on the inflation indicator and close stopcock.

2C



## POSITION BALLOON

- Grasp the device handle and align the device with the tissue tract.
- Pull gently to retract the device until the black line in the tension indicator window aligns with the markers on the side, indicating the balloon is abutting the venotomy with the correct amount of tension.

## DEPLOY THE SEALANT

3



## DEPLOY SEALANT

- While maintaining tension, press **button 1** until it is fully aligned with the handle, and the clock symbol is visible in the tension indicator window. This deploys and compresses the sealant against the venotomy.
- Lay the device down for **two minutes**.



## LOCK / STABILIZE / DEFLATE

4A



## LOCK SYRINGE

- Retract syringe plunger to lock position.




## STABILIZE VEIN

- Apply light fingertip compression proximal to the insertion site and then lightly grasp the device at skin with thumb and forefinger and realign with the tissue tract.

## DEFLATE BALLOON

- Open the stopcock to deflate the balloon.
- To ensure complete balloon deflation, wait until air bubbles and fluid have stopped moving through the inflation tubing.

## RETRACT THE BALLOON



## RETRACT BALLOON

- Pick up the device handle and realign with the tissue tract. Depress **button 2** to pull the deflated balloon into the device.

## REMOVE THE DEVICE

- While maintaining fingertip compression on the skin, remove the device from the patient.
- Continue to apply fingertip compression for up to **one minute** or as needed.
- Apply a sterile dressing once hemostasis is achieved.

# Ordering Information

**MYNX CONTROL™ VENOUS Vascular Closure Device** includes:

- (1) Balloon catheter with integrated sealant
- (1) 10 mL locking syringe

Size	Sheath Catch Color	Part Number
6F-12F	Black	MX61260

To order the MYNX CONTROL™ VENOUS Vascular Closure Device in the United States contact your local Cordis sales representative or customer service at 800.327.7714. To learn more visit

[cordis.com/mynxcontrolvenous](https://cordis.com/mynxcontrolvenous)

CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

IMPORTANT INFORMATION: Prior to use, refer to the Instructions for Use supplied with this device for indications, contraindications, suggested procedure, warnings, and precautions.

CORDIS, Cordis LOGO, and MYNX CONTROL are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2024 Cordis. All Rights Reserved. 100645248 06/2024.

