ENGLISH

SLEEK® OTW

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
ENGLISH

DESCRIPTION:

- The **SLEEK® OTW** Percutaneous Transluminal Angioplasty (PTA) Peripheral Catheter Family are a non-reusable coaxial design catheter with a semi-compliant balloon mounted on its distal tip.
- The hub/"Y" connector consists of a through lumen, allowing the catheter to track over a guidewire, and a balloon port, used to inflate the balloon (See diagram below).

- The catheter shaft is small beneath the balloon to achieve low profile. The distal tip is further tapered to accept the appropriate guidewire.


Figure 1


Figure 2

- Platinum iridium marker bands allow the balloon to be located under fluoroscopy. There are two marker bands denoting the ends of the working surface (Figure 1) on all sizes except the 1.25mm by 15mm and 1.5mm by 15mm balloons, where a single marker band is centered on the working surface (Figure 2).
- The SILX™ coating is a silicone coating which provides improved lubricity to the balloon shaft.
- A 0.014" (0.356 mm) guidewire is recommended for use with the **SLEEK® OTW** catheter.

Balloon Characteristics

- Please check the package label for the rated burst pressure (RBP). It is important that the balloon not be inflated beyond the rated burst pressure. Pressures in excess of rated burst pressure may cause the balloon to burst.

Compliance charts are provided with each product. Please note that balloon diameters may vary within manufacturing tolerances. All inflations should be viewed under fluoroscopy. The **SLEEK® OTW** balloons reach their nominal diameter at 6 atm (608 kPa).

INDICATIONS:

The **SLEEK® OTW** catheters are intended for balloon dilatation of the femoral, popliteal and infra-popliteal arteries. **These catheters are not designed to be used in the coronary arteries.**

WARNINGS:

- Contents supplied STERILE using ethylene oxide (EO). Non-pyrogenic. Do not reuse, reprocess or re-sterilize.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

- Visually inspect the packaging to verify that the sterile barrier is intact. Do not use if the sterile barrier is opened or damaged.
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- Do not exceed the rated burst pressure. A syringe with pressure gauge is recommended to monitor pressure. Pressure in excess of the rated burst pressure can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.
- Use a 20 ml or larger syringe for inflation.
- Use the catheter prior to the "Use By" date specified on the package.
- Do not advance the guidewire, balloon dilation catheter, or any component if resistance is met, without first determining the cause and taking remedial action.
- This catheter is not recommended for pressure measurement or fluid injection.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

PRECAUTIONS:

- Dilation procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Carefully inspect the catheter prior to use to verify that the catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
- Careful attention must be paid to the maintenance of tight catheter connections to avoid the introduction of air into the system.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath/guide catheter should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath/guide catheter as a unit and withdrawing both together, using a gently twisting motion combined with traction.
- Before removing catheter from sheath/guide catheter it is very important that the balloon is completely deflated.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.
- Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

ADVERSE EFFECTS:

Possible adverse effects include, but are not limited to, the following:

- Vessel perforation
- Vessel spasm
- Haemorrhage
- Haematoma
- Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection
- Systemic embolization
- Endocarditis

- Short-term hemodynamic deterioration
- Death
- Vascular thrombosis
- Drug reactions, allergic reaction to contrast media.
- Arteriovenous fistula
- Vessel dissection

STORAGE:

Store in a cool, dark, dry place. Use the catheter prior to the "Use By" date specified on the package.

DIRECTIONS FOR USE:

Inspection and Preparation

- Remove the balloon sleeve by first withdrawing the shipping mandrel slightly and then slowly removing the sleeve while holding the catheter as close to the balloon as possible.
- If any resistance is felt, or if any stretching of the catheter is observed while removing the balloon sleeve, the product should not be used.
- The catheter should then be inspected for bends, kinks or stretched portions. Do not use if product damage is evident.
- Prepare a mixture of contrast medium and normal saline as per normal procedure. (Recommended 25%/75%)
- Attach a stopcock and a 20 ml syringe half filled with the contrast solution to the balloon port.
- Point the syringe nozzle downward and aspirate until all air is removed from the balloon.
- Turn the stopcock off and maintain the vacuum in the balloon.
- Purge the catheter through lumen thoroughly.
- Reinserting the balloon into the balloon sleeve may damage the balloon or catheter.

PROCEDURE:

Insertion and Inflation Procedure:

NOTE: A 0.014" (0.356 mm) guidewire must be inserted in the SLEEK® OTW catheter across the balloon during any inflation of the balloon.

Make sure that the balloon sleeve has been removed from the catheter balloon. Enter the vessel percutaneously using the standard Seldinger technique over the appropriate guidewire for the size catheter being used.

Advance the catheter across the lesion with fluoroscopic guidance using accepted percutaneous transluminal angioplasty technique and inflate the balloon to the appropriate pressure.

NOTE: Do not inflate the balloon or advance the catheter unless the guidewire is in place.







Deflation and Withdrawal

- Deflate the balloon by drawing a vacuum with a 20 ml or larger syringe.
NOTE: The larger the syringe diameter, the greater the suction that is applied. For maximum deflation a 50 ml syringe is recommended.
- Gently withdraw the catheter. As the balloon exits the vessel, use a smooth, gentle, steady, counterclockwise motion. If resistance is felt upon removal then the balloon and the sheath should be removed together as a unit under fluoroscopic guidance, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.

Warning: After use this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practices and with applicable local, state and federal laws and regulations.

PRODUCT WARRANTY

We warrant that at the time of manufacture, these products were prepared and tested to verify that they are true to labelled claims. Due to biological differences in individuals, no product is 100% effective under all circumstances. In addition, because we have no control of the condition under which these products are used, diagnosis of the patient, the method of use or administration, and handling of the products after they leave our possession, we do not warrant for either a good effect or against ill effect following/during the use of the product. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES EITHER WRITTEN, ORAL OR IMPLIED (INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PURPOSE). No representative of the company may change any of the foregoing and the buyer accepts the product subject to all terms hereof.

EXPLANATION OF SYMBOLS ON LABELS AND PACKAGING			
	USE-BY DATE		LOT NUMBER
	STERILIZED USING ETHYLENE OXIDE		DO NOT REUSE
	NON-PYROGENIC		CONSULT INSTRUCTIONS FOR USE
	N UNITS PER BOX		MANUFACTURER
	OVER THE WIRE		KEEP DRY
	BALLOON PRESSURE		NOMINAL PRESSURE
	RECOMMENDED CATHETER SHEATH INTRODUCER		CAUTION
	BALLOON LENGTH		INFLATED BALLOON LENGTH IN mm
	BALLOON OUTER DIAMETER		INFLATED BALLOON DIAMETER IN mm
	USABLE CATHETER LENGTH		USABLE CATHETER LENGTH IN cm
	RECOMMENDED GUIDEWIRE		CATHETER DIAMETER

© 2017 Cardinal Health. All Rights Reserved. CORDIS, the Cordis LOGO, CARDINAL HEALTH and Sleek are trademarks of Cardinal Health and may be registered in the US and /or in other countries.

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Technologies Ltd