


Cordis[®]

EXOSEAL[®]


Vascular Closure Device

Precise and secure extravascular closure


Leave nothing behind¹




Precise & secure
closure



Trusted
bio-absorbable
technology




Easy to use



EXOSEAL[®]
Vascular Closure Device

Precise and secure extravascular closure



Deployment button
Lockout feature reduces risk of intravascular plug deployment

Sheath adaptor & wire loop release

- Compatible with introducer sheath (sheath exchange not required)²
- Minimises chance of infection and inflammation
- Maximises patient comfort

Indicator window
Provides indication of extravascular position before deployment

EXOSEAL[®]

EXOSEAL[®] Vascular Closure Device

French Size	Order Number	Description
5	EX500CE	EXOSEAL [®] - to be used with a 5F sheath introducer up to 12cm
6	EX600CE	EXOSEAL [®] - to be used with a 6F sheath introducer up to 12cm
7	EX700CE	EXOSEAL [®] - to be used with a 7F sheath introducer up to 12cm

REFERENCES: 1. Fully resorbed in 60-90 days. 2. Compatible with sheaths up to 12cm
IMPORTANT INFORMATION: Prior to use, refer to the Instructions for Use supplied with the device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous development, we reserve the right to change product specifications without prior notification.