

**Cordis**<sup>®</sup>

**MYNXGRIP**<sup>®</sup>


Vascular Closure Device

Secure extravascular closure without the trade offs<sup>1</sup>


Leave nothing behind<sup>1</sup>



Water soluble,  
bio-inert,  
non-thrombogenic  
polymer sealant



Secure adhesion to  
the vessel wall for  
active closure



Sealant swells to 3-4  
times original size filling  
the tissue tract



GRIP actively adheres  
to artery or vein providing  
durable haemostasis  
and a platform for  
natural healing




**MYNXGRIP®**

Vascular Closure Device

Secure extravascular closure  
without the trade offs<sup>1</sup>

Leave nothing behind<sup>1</sup>



MYNXGRIP® Vascular Closure Device

French Size	Colour	Order Number
5	Grey	MX5021
6/7	Green	MX6721

**Includes:** 1x Balloon catheter with integrated sealant, 1x 10ml locking syringe

**REFERENCES:** 1. The sealant is resorbed by the body within 30 days.

**IMPORTANT INFORMATION:** Prior to use, refer to the Instructions for Use supplied with the device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous development, we reserve the right to change product specifications without prior notification.