



# SABER<sup>®</sup>

PTA DILATATION CATHETER

セイバー

More .018" OTW Balloon Sizes  
Designed With More Crossability  
More Power. More Control.





カタログ番号	バルーン径 (mm)	バルーン有効長 (mm)	カテーテル有効長 (cm)	RBP (atm)	適合最小シースサイズ (F)	シャフト径 (F)	バルーンマーカー
48002002S	2.0	20	90	18	4	3.9	2
48002004S	2.0	40	90	18	4	3.9	2
48002010S *	2.0	100	90	18	4	3.9	3
48002015S *	2.0	150	90	18	4	3.9	3
48002504S *	2.5	40	90	18	4	3.9	2
48002510S *	2.5	100	90	18	4	3.9	3
48003002S	3.0	20	90	18	4	3.9	2
48003004S	3.0	40	90	18	4	3.9	2
48003010S *	3.0	100	90	18	4	3.9	3
48003015S *	3.0	150	90	14	4	3.9	3
48004002S	4.0	20	90	18	4	3.9	2
48004004S	4.0	40	90	18	4	3.9	2
48004010S	4.0	100	90	18	4	3.9	3
48004015S	4.0	150	90	14	4	3.9	3
48005002S	5.0	20	90	16	4	3.9	2
48005004S	5.0	40	90	16	4	3.9	2
48005010S	5.0	100	90	16	4	3.9	3
48005015S	5.0	150	90	14	5	3.9	3
48006002S	6.0	20	90	16	4	3.9	2
48006004S	6.0	40	90	16	4	3.9	2
48006010S	6.0	100	90	16	5	3.9	3
48006015S	6.0	150	90	14	5	3.9	3
48007002S	7.0	20	90	14	4	3.9	2
48007004S	7.0	40	90	14	4	3.9	2
48008002S	8.0	20	90	12	4	3.9	2
48008004S	8.0	40	90	12	4	3.9	2
48002002X	2.0	20	150	18	4	3.9	2
48002004X	2.0	40	150	18	4	3.9	2
48002010X *	2.0	100	150	18	4	3.9	3
48002015X *	2.0	150	150	18	4	3.9	3
48003002X	3.0	20	150	18	4	3.9	2
48003004X	3.0	40	150	18	4	3.9	2
48003010X *	3.0	100	150	18	4	3.9	3
48003015X *	3.0	150	150	14	4	3.9	3
48004002X	4.0	20	150	18	4	3.9	2
48004004X	4.0	40	150	18	4	3.9	2
48004010X	4.0	100	150	18	4	3.9	3
48004015X	4.0	150	150	14	4	3.9	3
48005002X	5.0	20	150	16	4	3.9	2
48005004X	5.0	40	150	16	4	3.9	2
48005010X	5.0	100	150	16	4	3.9	3
48005015X	5.0	150	150	14	5	3.9	3
48006002X	6.0	20	150	16	4	3.9	2
48006004X	6.0	40	150	16	4	3.9	2
48006010X	6.0	100	150	16	5	3.9	3
48006015X	6.0	150	150	14	5	3.9	3
48007002X *	7.0	20	150	14	4	3.9	2
48007004X *	7.0	40	150	14	4	3.9	2
48008002X *	8.0	20	150	12	4	3.9	2
48008004X *	8.0	40	150	12	4	3.9	2

\*印のついている製品は在庫をご確認の上、ご注文ください。

◆適合最大ガイドワイヤ径:0.018"

■販売名: SABER PTAカテーテル ■承認番号: 22600BZX00271000 ■クラス分類: IV ■JMDNコード: 17184014

