




SUPERTORQUE® PLUS
Diagnostic Catheters

Shapes you need,
quality you expect


Excellent pushability without the risk of kinking

Clear visibility with radiopaque tip




SUPERTORQUE® PLUS Diagnostic Catheters


Shapes you need, quality you expect



Body: Wire-braided polyurethane for precise torque control and maximum pushability



Distal Tip: Soft, radiopaque polyurethane tip for atraumatic, safe cannulation



Contact your local Cordis® Sales Representative or refer to the Catalogue for codes.

IMPORTANT INFORMATION: Prior to use, refer to the Instructions for Use supplied with the device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous development, we reserve the right to change product specifications without prior notification.

