




GET MORE

Control
Back-up support
Space to work



ADROIT™
GUIDING CATHETER






Hybrid design features two different sizes of wire braid

BACK-UP SUPPORT

- Engineered with Hybrid braided wire technology
- Controlled Back-up Support

Back-up Support^{1,2}



¹ Cordis 2012 and 2015 Data on file.


² Sample size was ≤ 3 .

CONTROL


Gain greater positioning confidence through:

- Hybrid braided wire technology that delivers excellent torque response
- Multi-segment design that provides kink resistance and minimizes whip

Torque Response^{1,2}



Whip^{1,2}




¹ Cordis 2012 and 2015 Data on file.

² Sample size was ≤ 3 .


SPACE TO WORK

Experience more room to maneuver with the available Larger Lumen ID (6F)*:

- .072" ID provides device compatibility and provides good visualization
- Large lumen to perform kissing balloon technique
- PTFE lining provides lubricious lumen for smooth delivery
- True lumen Technology ensures consistent ID from hub to tip




Kissing Balloon Insertion Force^{1,2}



Force to track second balloon (g)

Cordis ADROIT™ Guiding Catheter

Medtronic Launcher™ Guiding Catheter



Hybrid braided wire technology enables large lumen diameter (.072")

* Compared to Cordis Vista Brite Tip™ Guiding Catheter.


¹ Cordis 2012 and 2015 Data on file.

² For 6F configurations, sample size was ≤ 3.

DESIGNED FOR ATRAUMATIC PLACEMENT


Engineered to avoid vessel wall damage via:

- Longer, soft radiopaque tip provides enhanced visibility* and allows active cannulation.
- A multi-segment design that provides kink resistance and minimizes whip



* Compared to the Cordis VISTA BRITE TIP™ Guiding Catheter.

¹ The primary differences between the ADROIT™ Guiding Catheter and the VISTA BRITE TIP™ Guiding Catheter are the larger inner diameter of the ADROIT™ Guiding Catheter and the three transition segments in the ADROIT™ Guiding Catheter as opposed to two in the VISTA BRITE TIP™ Guiding Catheter.



ADROIT™ Guiding Catheter

Get More.

Control

A multi-segment design that provides kink resistance and minimizes whip

Back-up support

Hybrid braid construction delivers support, torque control and lumen integrity



Space to work

Large lumen to perform kissing balloon technique and PTFE lining provides lubricious lumen for smooth delivery


Shapes and Ordering Information


Cordis **ADROIT™** Guiding Catheter


A large portfolio of shapes and configurations are available to meet varying procedural and anatomical demands.


			PRODUCT CODE
Left Coronary Artery	Shape	Description	ADROIT™ 6F
	JL 3		672-000-00
	JL 3.5		672-002-00
	JL 4		672-004-00
	AL .75		672-034-00
	AL 1		672-036-00
	AL 2 LBT	Long BRITE TIP™	672-040-0L

Shapes and Ordering Information



Left Coronary Artery (cont.)	Shape	Description	PRODUCT CODE
			ADROiT™ 6F
<p>Extra Back-up</p> 	XB 3		672-052-00
	XB 3.25		672-400-00
	XB 3.5		672-054-00
	XB 3.5 SH	2 Side Holes	672-055-00
	XB 3.5	ECOPAK (5 units)	672-054-0E
	XB 4		672-056-00


<p>Extra Back-up LAD</p> 	XBLAD 3.5		672-060-00
	XBLAD 4		672-062-00

<p>Extra Back-up Contralateral</p> 	XBC 3.5		672-072-00
--	---------	--	------------


Right Coronary Artery	Shape	Description	PRODUCT CODE
			ADROiT™ 6F
<p>Judkins Right</p> 	JR 3.5		672-080-00
	JR 4		672-082-00
	JR 4 (90 cm)		672-082-90
	JR 4 (125 cm)		672-082-25
	JR 4 SH	2 Side Holes	672-083-00
	JR 4 LBT	Long BRITE TIP™	672-082-0L
	JR 4	ECOPAK (5 units)	672-082-0E
	JR 5		672-084-00

Shapes and Ordering Information


			PRODUCT CODE
Right Coronary Artery (cont.)	Shape	Description	ADRCiT™ 6F
	AR 1		672-110-00
	AR 1 SH	2 Side Holes	672-111-00
	AR 2		672-112-00
Extra Back-up RCA		XB RCA	672-126-00
3 DRC-Williams Technique		3 DRC	672-130-00

			PRODUCT CODE
Dedicated Transradial Shapes	Shape	Description	ADRCiT™ 6F
Right Radial Extra Back-up		RRAD6 XB	672-301-00

Shapes and Ordering Information

			PRODUCT CODE
Multipurpose	Shape	Description	ADROIT™ 6F
Multipurpose		MPA 1	672-270-00
		MPA 1 (125 cm)	672-270-25

Hockey Stick



H-STICK

672-278-00

			PRODUCT CODE
Bypass	Shape	Description	ADROIT™ 6F


Left Coronary Bypass



LCB

672-180-00


Right Coronary Bypass



RCB

672-182-00

Internal Mammary




IM

672-190-00

			PRODUCT CODE
Renal Curve	Shape	Description	ADROIT™ 6F

Renal Curve



RDC

672-212-55

ADROIT™ Guiding Catheter from Cordis

Speak with your Cordis representative about getting more out of your guiding catheter.

