

AVANTI™ + Sheath Introducer

For atraumatic insertion

Smooth performance
from start to finish




Cordis™

AVANTI™ + Sheath Introducer

For atraumatic insertion

The proven cannula offers reliable, smooth performance from start to finish



- For smooth, atraumatic insertion
- Kink resistant
- Optimized catheter maneuverability with minimal leakage

For smooth, atraumatic insertion

- Smooth transitions between cannula, guidewire and dilator for:
 - Effortless vessel entry
 - Increased patient comfort
 - Smallest possible risk for complications

Kink Resistant

Dual layer material gives AVANTI™ + Sheath Introducer excellent flexibility and the strength needed to prevent kinking - even in tortuous vessels and scar tissue

Optimized catheter maneuverability with minimal leakage

The patented SLIX haemostatic valve with six-cut spiral design allows free catheter movement without blood reflux or air aspiration



AVANTI™ + Sheath Introducer


For atraumatic insertion

Smooth Performance from Start to Finish

Available with or without mini-guidewire
(mid-length without wire only)

Available in three lengths
(usable cannula length)

- 5.5cm **brachial**
- 7.5cm **short length**
- 11cm **standard length**
- 23cm **mid-length**



CARDIOLOGY: AVANTI™ + Sheath Introducer

French size	Colour code	5.5cm Cannula		11cm Cannula		23cm Cannula
		With mini-wire	Without mini-wire	With mini-wire	Without mini-wire	Without mini-wire
4	■	402-604P		402-604X	402-604A	
5	■	402-605P		402-605X	402-605A	402-605T
5.5**	■			402-655X		
6	■	402-606P	402-606R	402-606X	402-606A	402-606T
6.5**	■			402-656X		
7	■	402-607P	402-607R	402-607X	402-607A	
8	■			402-608X	402-608A	402-608T
9	■			402-609X	402-609A	402-609T
10	■			402-610X	402-610A	
11	■			402-611X		

**

5F MS = Mid-Size = 5.5F
6F MS = Mid-Size = 6.5F

AVANTI™ + Sheath Introducer

For atraumatic insertion

CARDIOLOGY: AVANTI™+ Sheath Introducer

11 cm Cannula				
French size	Colour code	Mini-wire (inch)	With mini-wire	Without mini-wire
4	■	.035	504-604X	504-604A
5	■	.038	504-605X	504-605A
5.5**	■	.038	504-655X*	
6	■	.038	504-606X	504-606A
6.5**	■	.038	504-656X*	504-656A*
7	■	.038	504-607X	504-607A
7.5**	■	.038	504-657X*	
8	■	.038	504-608X	504-608A
8.5**	■	.038	504-658X*	
9	■	.038	504-609X	504-609A
10	■	.038	504-610X	504-610A
11	■	.038	504-611X	

* Monitoring sheath: provides slightly larger ID for pressure monitoring or side port aspiration with a catheter in place; they are one-half French size larger in size

** 5F MS = Mid-Size = 5.5F
 6F MS = Mid-Size = 6.5F
 7F MS = Mid-Size = 7.5F
 8F MS = Mid-Size = 8.5F



CARDIOLOGY: AVANTI™+ Sheath Introducer

23 cm Cannula				
French size	Colour code	Dilator usable (cm)	With obturator	Without obturator
4	■	30		504-604T
5	■	30		504-605T
6	■	30	504-606D	504-606T
6.5**	■	30		504-656T*
7	■	30	504-607D	504-607T
8	■	30	504-608D	504-608T
10	■	30		504-610T
11	■	30		504-611T

** 6F MS = Mid-Size = 6.5F



AVANTI™ + Sheath Introducer

For atraumatic insertion


CARDIOLOGY: AVANTI™ + Brachial Sheath Introducer with Mini-Guidewire

5.5 cm Cannula				
French size	Colour code	Mini-wire (inch)	Catalogue number	
4	■	.035	504-604P	
5	■	.038	504-605P	
6	■	.038	504-606P	



CARDIOLOGY: AVANTI™ + Short length Sheath Introducer with mini-guidewire

7.5 cm Cannula				
French size	Colour code	Mini-wire (inch)	Catalogue number	
4	■	.021	504-604S	
5	■	.021	504-605S	
6	■	.021	504-606S	



CARDIOLOGY: AVANTI™ + Transradial Sheath Introducer kit with Mini-Guidewire

French size	Colour code	11 cm Cannula		23 cm Cannula	
		Wire acceptance (inch)	Catalogue number	Wire acceptance (inch)	Catalogue number
4	■	.021	504-614Z	.021 70 cm	504-624Z
5	■	.021	504-615Z	.021 70 cm	504-625Z
6	■	.021	504-616Z	.021 70 cm	504-626Z
7	■	.021	504-617Z		

