


S.M.A.R.T.™ Flex

Vascular Stent System

**Our inspiration came
right from the source.**



Experience our Cordis
self-expanding stents.

S.M.A.R.T.™ Flex Vascular Stent System

Feeling is believing

Cordis™

S.M.A.R.T.™ Flex Vascular Stent System Product Description

Type: Fully connected yet flexible stent design
Material: Nitinol, with tantalum markers on stent ends
Maximum Guidewire: 0.035"
Stent Lengths: 20-200 mm
Sheath Compatibility: 6F

Stent Diameters: 5-10 mm
Stent Delivery System Working Lengths: 80 cm (S suffix) and 120 cm (M suffix)

Ordering Information

Stent Diameter (mm)	Stent Length (mm)	Recommended Lumen Size (mm)	Recommended Sheath Introducer Size	Catalogue Number	
				Length of delivery system 80 cm	Length of delivery system 120 cm
5	30	3.5-4.5	6	SF05030SV	SF05030MV
5	40	3.5-4.5	6	SF05040SV	SF05040MV
5	60	3.5-4.5	6	SF05060SV	SF05060MV
5	80	3.5-4.5	6	SF05080SV	SF05080MV
5	100	3.5-4.5	6	SF05100SV	SF05100MV
5	120	3.5-4.5	6	SF05120SV	SF05120MV
5	150	3.5-4.5	6	SF05150SV	SF05150MV
6	30	4.5-5.5	6	SF06030SV	SF06030MV
6	40	4.5-5.5	6	SF06040SV	SF06040MV
6	60	4.5-5.5	6	SF06060SV	SF06060MV
6	80	4.5-5.5	6	SF06080SV	SF06080MV
6	100	4.5-5.5	6	SF06100SV	SF06100MV
6	120	4.5-5.5	6	SF06120SV	SF06120MV
6	150	4.5-5.5	6	SF06150SV	SF06150MV
7	30	5.5-6.5	6	SF07030SV	SF07030MV
7	40	5.5-6.5	6	SF07040SV	SF07040MV
7	60	5.5-6.5	6	SF07060SV	SF07060MV
7	80	5.5-6.5	6	SF07080SV	SF07080MV
7	100	5.5-6.5	6	SF07100SV	SF07100MV
7	120	5.5-6.5	6	SF07120SV	SF07120MV
7	150	5.5-6.5	6	SF07150SV	SF07150MV
7	200	5.5-6.5	6	SF07200SV	SF07200MV
8	30	6.5-7.5	6	SF08030SV	SF08030MV
8	40	6.5-7.5	6	SF08040SV	SF08040MV
8	60	6.5-7.5	6	SF08060SV	SF08060MV
8	80	6.5-7.5	6	SF08080SV	SF08080MV
8	100	6.5-7.5	6	SF08100SV	SF08100MV
8	120	6.5-7.5	6	SF08120SV	SF08120MV
8	150	6.5-7.5	6	SF08150SV	SF08150MV
8	200	6.5-7.5	6	SF08200SV	SF08200MV



For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and S.M.A.R.T. are trademarks of Cordis and may be registered in the US and/or in other countries.
 © 2023 Cordis. All Rights Reserved. 100551514-4 01/2023