


# Cordis EMERALD™ Diagnostic Guidewires



PTFE pre-coated  
guidewire

Rosen Heavy Duty and Amplatz  
Super Stiff EMERALD™ Guidewires

# Cordis EMERALD™

## Guidewires

### Rosen Heavy Duty and Amplatz Super Stiff PTFE-Coated EMERALD™ Guidewires\*

#### EMERALD™ Amplatz Super Stiff Guidewires

Diameter (inches)	Length (cm)	J-Radius (mm)	Catalogue Number
PTFE-Coated STRAIGHT TIP			
.035 Amplatz	150	–	502-726
.035 Amplatz	180	–	502-728

#### PTFE-Coated J-CURVE

.035 Amplatz	150	3	502-731
.035 Amplatz	180	3	502-733
.035 Amplatz	260	3	502-735

- For percutaneous entry and guidance of angiographic catheters
- Pre-coated PTFE
- 5 units per package

#### EMERALD™ Rosen Heavy Duty Guidewires

Diameter (inches)	Length (cm)	J-Radius (mm)	Catalogue Number
PTFE-Coated J-CURVE			
.035 Rosen	180	1.5	502-717

- Atraumatic end
- For percutaneous entry and high support guidance of angiographic catheters
- 5 units per package

\* More configurations available. Please consult your Cordis Catalogue.

