

DESIGNED FOR EFFICIENCY







SABER™ .035
PTA DILATATION CATHETER




EXCELLENT NAVIGATION THROUGH TORTUOUS ANATOMY AND DIFFICULT LESIONS

KEY FEATURES FOR ENHANCED PERFORMANCE

 <p>Smooth tapered tip enhances crossability</p>	 <p>Dual layer hydrophilic coating improves trackability</p>	 <p>Dual distal markers on ≥ 100 mm balloons enhance visualization under fluoroscopy</p>	 <p>Proprietary balloon design offers excellent rewrap and facilitates smooth retraction through the sheath</p>
---	---	--	--

BALLOON CROSSABILITY RESULTS¹



Crossability testing performed with 6 x 200 mm balloons in a Superficial Femoral Artery (SFA) Crossing Model

OPTIMAL DESIGN FOR STREAMLINED CARE IN MORE PATIENTS

EXTENSIVE SIZE MATRIX EXPANDS TREATMENT APPLICATIONS

SABER™ .035 CATHETER SIZE MATRIX CHART* USABLE CATHETER LENGTHS: 80, 135 CM

		BALLOON LENGTH (MM)									
		20	40	60	80	100	120	150	200	250	300
BALLOON DIAMETER (MM)	3	24	24	24	22	22	22	22	22	20	20
	4	22	22	22	22	22	22	22	22	20	20
	5	22	22	22	22	22	20	20	20	18	18
	6	18	18	17	17	17	17	17	17	17	17
	7	18	18	16	16	16	14	14	14	14	14
	8	16	16	16	16	16					
	9	13	13	13	13	13					
	10	14	14	14	14	14					
	12	12	12								

Values in table represent RBP

5F	6F	7F
----	----	----

*Not all balloon sizes are available in all shaft lengths. Please consult the product catalog for additional information.

PORTFOLIO THAT ENHANCES RESULTS

<p>Balloon Lengths of up to 300 mm</p> <p>enable treatment of long lesions with fewer inflations and as a result may potentially reduce procedural and fluoroscopy time</p>	<p>Low Profile 5F Sheath Compatibility</p> <p>for most sizes may reduce the potential for access site complications</p>	<p>High Rated Burst Pressure</p> <p>of up to 24 atm for effective dilation of tough lesions</p>	<p>Comprehensive Product Portfolio</p> <p>reduces the need for multiple vendors</p>
--	--	--	--

SABER™.035

PTA DILATATION CATHETER

OPTIMIZES OUTCOMES WITH:



Superior Crossability¹



Excellent Deliverability



Broad Size Matrix



Contact your Cordis representative to evaluate SABER™.035 Catheters or for more information on the full line of lower-extremity solutions, visit www.cordis.com/emea/home

Follow & share YOUR experience on X & LinkedIn with #SABER35

1. Cordis 2020 data on file

