

NAVIGATE THE MOST
CHALLENGING
PATHWAYS WITH
LUBRICITY AND
CONTROL



TEMPO AQUA™
HYDROPHILIC-COATED DIAGNOSTIC CATHETER



Cordis TEMPO AQUA™ Diagnostic Catheter

The **TEMPO AQUA™** Catheter makes endovascular procedures quick and easy by combining the Cordis TEMPO™ Angiographic Catheter characteristics with a durable hydrophilic coating for easy advancement into the vasculature and subintimal space.
















CROSSABILITY














- Hydrophilic coating to **navigate the most challenging cases** with lubricity and control
- Hydrophilic coated distal shaft enables **smooth passage through the vessels**
- Catheter tip demonstrates superb shape retention, greater positional stability and is made from a special blend of soft **atraumatic** material

CONTROL

- Distal half of the catheter is coated for high lubricity, **proximal part is uncoated** for secure handling
- High-strength nylon construction and braided shaft design give strength and flexibility for **superb torque response** and maximum **tactile control**
- Thin-wall technology guarantees a large inner diameter and provides **high flow rates at low injection force**

ORDERING INFORMATION

TEMPO AQUA™ Catheter 4F				
Catalog Number	Description	Shape	Side Holes	Length, cm
SRD6913	 Straight	STR	0	65
SRD6821	 Straight	STR	0	100
SRD6667	 Straight	STR	0	125
452-407H0	 Multipurpose Adult	MPA	0	100
452-442V0	 Child 1	C1	0	65
452-443V0	 Small Adult	C2	0	65
452-443S0	 Small Adult	C2	0	80
452-430H0	 Sidewinder Simmons Technique I	SIM1	0	100
452-431H0	 Sidewinder Simmons Technique II	SIM2	0	100
452-435H0	 Headhunter I	H1	0	100
SRD6107	 Headhunter I	H1	0	125
452-413H0	 Berenstein	BER	0	100
SRD6912	 Vertebral	VERT	0	65
452-414H0	 Vertebral	VERT	0	100
SRD6108	 Vertebral	VERT	0	125

TEMPO AQUA™ Catheter 5F				
Catalog Number	Description	Shape	Side Holes	Length, cm
SRD6109	 Sidewinder Simmons Technique II	SIM2	0	125
SRD6669	 Straight	STR	0	125
452-507H0	 Multipurpose Adult	MPA	0	100
452-542V0	 Child 1	C1	0	65
452-543V0	 Small Adult	C2	0	65
SRD6677	 Small Adult	C2	0	80
452-530H0	 Sidewinder Simmons Technique I	SIM1	0	100
452-531H0	 Sidewinder Simmons Technique II	SIM2	0	100
452-535H0	 Headhunter I	H1	0	100
SRD6721	 Headhunter I	H1	0	125
452-513H0	 Berenstein	BER	0	100
SRD6906	 Vertebral	VERT	0	65
452-514H0	 Vertebral	VERT	0	100

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO, TEMPO and TEMPO AQUA are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2024 Cordis. All Rights Reserved. 100652010 09/2024

CE
2797

