
BIOPSY FORCEPS



Cordis.

Cordis Biopsy Forceps

EXCELLENT QUALITY SAMPLES

- Two stainless steel symmetrical hinged cutting jaws.
- Myocardial tissues cut cleanly for better histological analysis.

ADAPTABILITY TO A WIDE RANGE OF CLINICAL SITUATIONS




Jaw volume = 1.84mm^3

REDUCED RISK OF TISSUE TRAUMA DURING SAMPLING



Cordis Biopsy Forceps



ADAPTABILITY TO A WIDE RANGE OF CLINICAL SITUATIONS

Choice of two sizes of forceps diameter for taking samples adapted to a wide range of clinical situations.

- 2.46 mm³ of tissue sample using the 5.5F forceps
- 5.20 mm³ of tissue sample using the 7F forceps Jugular or femoral access possible with two lengths of forceps.



EXCELLENT QUALITY SAMPLES

- Two stainless steel symmetrical hinged cutting iaws.
- Myocardial tissues cut cleanly for better histological analysis.
- Reduced risk of tissue trauma



GREAT OPERATIONAL FLEXIBILITY

- 3-pull ring handle with spring so that samples can be held easily while being removed.
- Central ring of the handle rotates to accommodate any hand position.

BIOPSY FORCEPS - Ordering information

| Standard Biopsy Forceps | | | | | | |
|--------------------------|--------------------|---------------|-------------------------------|-------------------|--|--|
| Shaft OD, French Size | Shaft OD, mm | Length, cm | Description | Catalog Number | | |
| 5.5 | 1.85 | 104 | For femoral approach | 504-300 | | |
| 5.5 | 1.85 | 50 | For internal jugular approach | 504-302 | | |
| 7.0 | 2.3 | 104 | For femoral approach | 504-300L | | |
| 7.0 | 2.3 | 50 | For internal jugular approach | 504-302L | | |

| BIPAL™ Biopsy Forceps | | | | | | |
|--------------------------|--------------------|---------------|--|-------------------|--|--|
| Shaft OD, French Size | Shaft OD, mm | Length, cm | Description | Catalog Number | | |
| 7.0 | 2.3 | 104 | Straight tip for femoral approach | 502-400B | | |
| 7.0 | 2.3 | 50 | Radial tip shape for internal jugular approach | 502-402B | | |
| 7.0 | 2.3 | 50 | Multipurpose tip shape for internal jugular approach | 502-402M | | |

| Sheath Introducers and Sheath Sets | | | | | | |
|---|--------------------------|---------------------------------|-------------------|--|--|--|
| For use with Catheter French Size | Sheath, Length, cm | Description | Catalog Number | | | |
| 7 | 43 | Sheath only, Straight tip | 501-611 | | | |
| 7 | 96 | Sheath only, Straight tip | 501-613 | | | |
| 7 | 43 | Sheath only, Multipurpose curve | 501-611A | | | |
| 7 | 96 | Sheath only, Multipurpose curve | 501-613A | | | |

