


ADVANCING
THE CUTTING
EDGE IN BIOPSY
FORCEPS



BIOPSY FORCEPS



Cordis Biopsy Forceps



ADAPTABILITY TO A WIDE RANGE OF CLINICAL SITUATIONS

Choice of two sizes of forceps diameter for taking samples adapted to a wide range of clinical situations.

- 2.46 mm³ of tissue sample using the 5.5F forceps
- 5.20 mm³ of tissue sample using the 7F forceps

Jugular or femoral access possible with two lengths of forceps.



EXCELLENT QUALITY SAMPLES

- Two stainless steel symmetrical hinged cutting jaws.
- Myocardial tissues cut cleanly for better histological analysis.
- Reduced risk of tissue trauma



GREAT OPERATIONAL FLEXIBILITY

- 3-pull ring handle with spring so that samples can be held easily while being removed.
- Central ring of the handle rotates to accommodate any hand position.

BIOPSY FORCEPS - Ordering information

Standard Biopsy Forceps				
Shaft OD, French Size	Shaft OD, mm	Length, cm	Description	Catalog Number
5.5	1.85	104	For femoral approach	504-300
5.5	1.85	50	For internal jugular approach	504-302
7.0	2.3	104	For femoral approach	504-300L
7.0	2.3	50	For internal jugular approach	504-302L


BIPAL™ Biopsy Forceps				
Shaft OD, French Size	Shaft OD, mm	Length, cm	Description	Catalog Number
7.0	2.3	104	Straight tip for femoral approach	502-400B
7.0	2.3	50	Radial tip shape for internal jugular approach	502-402B
7.0	2.3	50	Multipurpose tip shape for internal jugular approach	502-402M

Sheath Introducers and Sheath Sets			
For use with Catheter French Size	Sheath Length, cm	Description	Catalog Number
7	43	Sheath only, Straight tip	501-611
7	96	Sheath only, Straight tip	501-613
7	43	Sheath only, Multipurpose curve	501-611A
7	96	Sheath only, Multipurpose curve	501-613A


502-400B



502-402B



502-402M



For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO and BIPAL are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2024 Cordis. All Rights Reserved. 100557686-6 09/2024

