


Cordis **MAXI LD**[™]

PTA Dilatation Catheter

The 0.035"
over-the-wire system
of choice for treating
peripheral arteries.



Cordis **MAXI LD™**

PTA Dilatation Catheter

MAXI LD™

The 0.035" over-the-wire system of choice for treating peripheral arteries.

Ordering Information

Balloon OD x Length (mm x cm)	Rated Burst Pressure (atm)	Recommended Sheath Fit (F)	Ordering Number	Usable Catheter Length (cm)	
				80	110
14 x 4	6	8	416-1440	S	L
14 x 6	6	8	416-1460	S	L
15 x 4	6	9	416-1540	S	L
15 x 6	6	9	416-1560	S	L
16 x 4	6	10	416-1640	S	L
16 x 6	6	10	416-1660	S	L
18 x 4	6	10	416-1840	S	L
18 x 6	6	10	416-1860	S	L
20 x 4	6	11	416-2040	S	L
20 x 6	6	11	416-2060	S	L
20 x 8	6	11	416-2080		LF
22 x 4	5	12	416-2240	S	L
25 x 4	5	12	416-2540	S	L

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and MAXI LD are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2020 Cardinal Health. All Rights Reserved. 100557571 03/2020

For more information please
visit www.cardinalhealth.co.uk



Cordis™
A Cardinal Health company