

# PALMAZ GENESIS™ Stent on OPTA™ Pro PTA Catheter

(Large)

## .035" Over-the-wire System

PALMAZ GENESIS™ Stent on OPTA™ Pro (Large) PTA Catheter

- Recommended inflation pressure : 8 atm

## Ordering Information

| Balloon<br>diameter x length<br>(mm x mm) | Unexpanded<br>Stent Length<br>(mm) | Expanded<br>Stent Length<br>(mm) | Sheath<br>Compatibility<br>(F) | Guiding<br>Compatibility<br>(F) | Catalog Number<br>Length of Delivery System |           |
|---|------------------------------------|----------------------------------|--------------------------------|---------------------------------|---|-----------|
|   |                                    |                                  |                                |                                 | 80 cm                                       | 135 cm    |
| 5 x 40                                    | 39                                 | 39                               | 6                              | 8 (.088")                       | PG3950PPS                                   | PG3950PPX |
| 5 x 60                                    | 59                                 | 59                               | 6                              | 8 (.088")                       | PG5950PPS                                   | PG5950PPX |
| 6 x 30                                    | 29                                 | 27                               | 6                              | 8 (.088")                       | PG2960PPS                                   | PG2960PPX |
| 8 x 30                                    | 29                                 | 26                               | 7                              | 9 (.098")                       |   | PG2980PPX |
| 6 x 40                                    | 39                                 | 38                               | 6                              | 8 (.088")                       | PG3960PPS                                   | PG3960PPX |
| 6 x 60                                    | 59                                 | 59                               | 6                              | 8 (.088")                       | PG5960PPS                                   | PG5960PPX |
| 6 x 80                                    | 79                                 | 79                               | 6                              | 8 (.088")                       | PG7960PPS                                   | PG7960PPX |
| 7 x 30                                    | 29                                 | 26                               | 6                              | 8 (.088")                       | PG2970PPS                                   | PG2970PPX |
| 7 x 40                                    | 39                                 | 36                               | 6                              | 8 (.088")                       |   | PG3970PPX |
| 7 x 60                                    | 59                                 | 57                               | 6.5                            | 8 (.088")                       | PG5970PPS                                   | PG5970PPX |
| 7 x 80                                    | 79                                 | 78                               | 6.5                            | 8 (.088")                       | PG7970PPS                                   | PG7970PPX |
| 8 x 40                                    | 39                                 | 36                               | 7                              | 9 (.098")                       | PG3980PPS                                   | PG3980PPX |
| 8 x 60                                    | 59                                 | 56                               | 7                              | 9 (.098")                       | PG5980PPS                                   | PG5980PPX |
| 8 x 80                                    | 79                                 | 76                               | 7                              | 9 (.098")                       | PG7980PPS                                   | PG7980PPX |
| 9 x 20                                    | 19                                 | 17                               | 7                              | 9 (.098")                       | PG1990PPS                                   | PG1990PPX |
| 9 x 26                                    | 25                                 | 22                               | 7                              | 9 (.098")                       | PG2590PPS                                   | PG2590PPX |
| 9 x 30                                    | 29                                 | 25                               | 7                              | 9 (.098")                       | PG2990PPS                                   | PG2990PPX |
| 9 x 40                                    | 39                                 | 35                               | 7                              | 9 (.098")                       | PG3990PPS                                   | PG3990PPX |
| 9 x 60                                    | 59                                 | 53                               | 7                              | 9 (.098")                       | PG5990PPS                                   | PG5990PPX |
| 9 x 80                                    | 79                                 | 75                               | 7                              | 9 (.098")                       | PG7990PPS                                   | PG7990PPX |
| 10 x 20                                   | 19                                 | 16                               | 7                              | 9 (.098")                       | PG1910PPS                                   | PG1910PPX |
| 10 x 26                                   | 25                                 | 22                               | 7                              | 9 (.098")                       |   | PG2510PPX |
| 10 x 30                                   | 29                                 | 24                               | 7                              | 9 (.098")                       | PG2910PPS                                   | PG2910PPX |
| 10 x 40                                   | 39                                 | 32                               | 7                              | 9 (.098")                       | PG3910PPS                                   | PG3910PPX |


# Complementary Cordis Products

The PALMAZ GENESIS™ Stent on OPTA™ Pro PTA Catheter can be used with other Cordis products. The products listed hereunder are only a selection. Please contact your local sales representative for further information on product compatibilities.

## Ordering Information

### .035" STORQ™

Steerable Guidewires



| Diameter (in) | Tip Flexibility | Tip Shape  | Length (cm) | Catalogue Number |
|---------------|-----------------|------------|-------------|------------------|
| .035          | Standard        | Straight   | 180         | 503-456          |
| .035          | Standard        | Angled     | 180         | 503-456J         |
| .035          | Standard        | Modified J | 180         | 503-456MJ        |
| .035          | Standard        | Straight   | 300         | 503-456X         |
| .035          | Standard        | Angled     | 300         | 503-456Y         |
| .035          | Standard        | Modified J | 300         | 503-456MY        |

### BRITE TIP™

Sheath Introducers



| Cannula usable length (cm) | Catalogue Number 6F | Catalogue Number 7F | Catalogue Number 8F |
|----------------------------|---------------------|---------------------|---------------------|
| 11                         | 401-611 M           | 401-711 M           | 401-811 M           |
| 23                         | 401-623 M           | 401-723 M           | 401-823 M           |
| 35                         | 401-635 M           | 401-735 M           | 401-835 M           |
| 45                         | 401-645 M           | 401-745 M           | 401-845 M           |
| 55                         | 401-655 M           | 401-755 M           | 401-855 M           |
| 90                         | 401-690 M           | 401-790 M           | 401-890 M           |