



THINK SMART

120 & 150
mm length



ALL LENGTHS COMPATIBLE WITH 6F SHEATH

S.M.A.R.T. CONTROL™
SELF EXPANDING NITINOL STENT





The longer S.M.A.R.T.™ Nitinol Stent System

Since their introduction in 1999, S.M.A.R.T.™ Stents have built and expanded their legacy. Every S.M.A.R.T.™ Stent is a reaffirmation of the Cordis Endovascular commitment to better patient care through evidence-based, results-oriented innovation.


120 & 150
mm length

WE'RE NOT JUST INCREASING THE LENGTH. WE'RE GROWING THE LEGACY.



A **legacy** of proven design

- **EXCELLENT SCAFFOLDING**
for uniform coverage
- **OFFSET ALTERNATING BRIDGES**
for wall apposition on a bend
- **SHORT STRUT LENGTH**
for individual cell response



A **legacy** of confidence

- **MORE THAN 1,000,000 STENT PLACEMENTS WORLDWIDE**
- **ROBUST CLINICAL EVIDENCE⁽¹⁻²⁻³⁾**
- **EXTENSIVELY STUDIED NITINOL STENT⁽¹⁻²⁻³⁾**

¹SIROCCO II, J. Endovasc. Ther. 2006; 13:701-710
²CRISP_US, J. Vasc. Interv. Radiol. 2004; 15:911-918
³STROLL- 3 year results, 2013

Think **SMART**

Cordis S.M.A.R.T.™ Nitinol Stent System

Long 120 mm and 150 mm Stent

ORDERING INFORMATION

Cordis **S.M.A.R.T.™** Nitinol Stent System
120 mm and 150 mm Stent Lengths

120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C06120MV	6	120	6/8
C07120MV	7	120	6/8
C08120MV	8	120	6/8
C06150MV	6	150	6/8
C07150MV	7	150	6/8

Cordis **S.M.A.R.T. CONTROL™** Nitinol Stent System
Large, 12 x 14 mm diameter

80cm Delivery System	120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C12040SV	C12040MV	12	40	7/9
C12060SV	C12060MV	12	60	7/9
C12080SV	C12080MV	12	80	7/9
C14040SV	C14040MV	14	40	7/9
C14060SV	C14060MV	14	60	7/9
C14080SV	C14080MV	14	80	7/9

Cordis **S.M.A.R.T. CONTROL™** Nitinol Stent System

80cm Delivery System	120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C06040SV	C06040MV	6	40	6/8
C06060SV	C06060MV	6	60	6/8
C06080SV	C06080MV	6	80	6/8
C06100SV	C06100MV	6	100	6/8
C07040SV	C07040MV	7	40	6/8
C07060SV	C07060MV	7	60	6/8
C07080SV	C07080MV	7	80	6/8
C07100SV	C07100MV	7	100	6/8
C08020SV		8	20	6/8
C08030SV		8	30	6/8
C08040SV	C08040MV	8	40	6/8
C08060SV	C08060MV	8	60	6/8
C08080SV	C08080MV	8	80	6/8
C08100SV	C08100MV	8	100	6/8
C09030SV		9	30	6/8
C09040SV	C09040MV	9	40	6/8
C09060SV	C09060MV	9	60	6/8
C09080SV	C09080MV	9	80	6/8
C10040SV	C10040MV	10	40	6/8
C10060SV	C10060MV	10	60	6/8
C10080SV	C10080MV	10	80	6/8

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO, S.M.A.R.T. CONTROL and S.M.A.R.T. are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2026 Cordis. All Rights Reserved. 100620931-2

