


### PRODUCT DESCRIPTION

**Type:** Low profile, high pressure, long shaft, rapid-exchange PTA balloon catheter designed to dilate stenoses in iliac, superficial femoral and proximal popliteal arteries via radial access. This device is also intended for post-dilatation of balloon-expandable and self-expanding stents.

<b>RBP RANGE</b>	<b>10-18 atm</b>	<b>SHEATH FIT RANGE</b>	<b>4F-6F</b>
<b>NOMINAL PRESSURE</b>	<b>6-8 atm</b>	<b>BALLOON MATERIAL</b>	<b>Nylon (Vestamid)</b>
<b>BALLOON LENGTHS</b>	<b>20-300 mm</b>	<b>MARKERS</b>	<b>Proximal and distal marker bands (dual distal marker bands for balloon lengths ≥100 mm)</b>
<b>BALLOON DIAMETERS</b>	<b>4-10 mm</b>	<b>GUIDEWIRE</b>	<b>.018"</b>
<b>SHAFT LENGTHS</b>	<b>190 cm</b>	<b>COATING</b>	<b>Dual layer hydrophilic coating</b>
<b>SHAFT MATERIAL</b>	<b>Pebax</b>		



### Ordering Information

DIAMETER (mm)	LENGTH (mm)	SHAFT LENGTH		
		RBP	SHEATH	190 cm
4	20	18	4	R51004002L
4	30	18	4	R51004003L
4	40	18	4	R51004004L
4	60	18	4	R51004006L
4	80	18	4	R51004008L
4	100	18	4	R51004010L
4	150	14	4	R51004015L
4	200	14	4	R51004020L
4	250	14	4	R51004025L
4	300	14	4	R51004030L

DIAMETER (mm)	LENGTH (mm)	SHAFT LENGTH		
		RBP	SHEATH	190 cm
5	20	16	4	R5100500L
5	30	16	4	R51005003L
5	40	16	4	R51005004L
5	60	16	4	R51005006L
5	80	16	4	R51005008L
5	100	16	4	R51005010L
5	150	14	4	R51005015L
5	200	14	4	R51005020L
5	250	14	4	R51005025L
5	300	14	4	R51005030L

## Ordering Information

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH	
					190 cm
6	20	16	4		R51006002L
6	30	16	4		R51006003L
6	40	16	4		R51006004L
6	60	16	4		R51006006L
6	80	16	5		R51006008L
6	100	16	5		R51006010L
6	150	14	5		R51006015L
6	200	14	5		R51006020L
6	250	14	6		R51006025L
6	300	14	6		R51006030L
7	20	14	4		R51007002L
7	30	14	4		R51007003L
7	40	14	5		R51007004L
7	60	14	5		R51007006L
7	80	14	5		R51007008L
7	100	14	5		R51007010L
8	20	12	4		R51008002L
8	30	12	4		R51008003L
8	40	12	4		R51008004L
8	60	12	5		R51008006L
8	80	12	5		R51008008L
8	100	12	5		R51008010L
9	20	10	5		R51009002L
9	30	10	5		R51009003L
9	40	10	5		R51009004L
9	60	10	5		R51009006L
9	80	10	5		R51009008L
9	100	10	5		R51009010L
10	20	10	5		R51010002L
10	30	10	5		R51010003L
10	40	10	5		R51010004L
10	60	10	5		R51010006L
10	80	10	6		R51010008L
10	100	10	6		R51010010L

For Healthcare Professionals Only.

**Important information:** Prior to use, refer to the instructions for use supplied with this device for indications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO, BRITE TIP RADIANTZ and SABERX RADIANTZ are trademarks of Cordis and may be registered in the US and/or in other countries. © 2025 Cordis. All Rights Reserved. 100677382

