


# S.M.A.R.T.™ Flex

## Vascular Stent System

**Our inspiration came  
right from the source.**



Experience our Cordis  
self-expanding stents.

**S.M.A.R.T.™ Flex Vascular Stent System**

**Feeling is believing**

**Cordis**

# S.M.A.R.T.™ Flex Vascular Stent System Product Description

**Type:** Fully connected yet flexible stent design

**Material:** Nitinol, with tantalum markers on stent ends

**Maximum Guidewire:** 0.035"

**Stent Lengths:** 20-200 mm

**Sheath Compatibility:** 6F

**Stent Diameters:** 5-10 mm

**Stent Delivery System Working Lengths:** 80 cm (S suffix) and 120 cm (M suffix)

## Ordering Information

| Stent Diameter (mm) | Stent Length (mm) | Recommended Lumen Size (mm) | Recommended Sheath Introducer Size | Catalogue Number<br>Length of delivery system |           |
|---------------------|-------------------|-----------------------------|------------------------------------|-----------------------------------------------|-----------|
|                     |                   |                             |                                    | 80 cm                                         | 120 cm    |
| 5                   | 30                | 3.5-4.5                     | 6                                  | SF05030SV                                     | SF05030MV |
| 5                   | 40                | 3.5-4.5                     | 6                                  | SF05040SV                                     | SF05040MV |
| 5                   | 60                | 3.5-4.5                     | 6                                  | SF05060SV                                     | SF05060MV |
| 5                   | 80                | 3.5-4.5                     | 6                                  | SF05080SV                                     | SF05080MV |
| 5                   | 100               | 3.5-4.5                     | 6                                  | SF05100SV                                     | SF05100MV |
| 5                   | 120               | 3.5-4.5                     | 6                                  | SF05120SV                                     | SF05120MV |
| 5                   | 150               | 3.5-4.5                     | 6                                  | SF05150SV                                     | SF05150MV |
| 6                   | 30                | 4.5-5.5                     | 6                                  | SF06030SV                                     | SF06030MV |
| 6                   | 40                | 4.5-5.5                     | 6                                  | SF06040SV                                     | SF06040MV |
| 6                   | 60                | 4.5-5.5                     | 6                                  | SF06060SV                                     | SF06060MV |
| 6                   | 80                | 4.5-5.5                     | 6                                  | SF06080SV                                     | SF06080MV |
| 6                   | 100               | 4.5-5.5                     | 6                                  | SF06100SV                                     | SF06100MV |
| 6                   | 120               | 4.5-5.5                     | 6                                  | SF06120SV                                     | SF06120MV |
| 6                   | 150               | 4.5-5.5                     | 6                                  | SF06150SV                                     | SF06150MV |
| 7                   | 30                | 5.5-6.5                     | 6                                  | SF07030SV                                     | SF07030MV |
| 7                   | 40                | 5.5-6.5                     | 6                                  | SF07040SV                                     | SF07040MV |
| 7                   | 60                | 5.5-6.5                     | 6                                  | SF07060SV                                     | SF07060MV |
| 7                   | 80                | 5.5-6.5                     | 6                                  | SF07080SV                                     | SF07080MV |
| 7                   | 100               | 5.5-6.5                     | 6                                  | SF07100SV                                     | SF07100MV |
| 7                   | 120               | 5.5-6.5                     | 6                                  | SF07120SV                                     | SF07120MV |
| 7                   | 150               | 5.5-6.5                     | 6                                  | SF07150SV                                     | SF07150MV |
| 7                   | 200               | 5.5-6.5                     | 6                                  | SF07200SV                                     | SF07200MV |
| 8                   | 30                | 6.5-7.5                     | 6                                  | SF08030SV                                     | SF08030MV |
| 8                   | 40                | 6.5-7.5                     | 6                                  | SF08040SV                                     | SF08040MV |
| 8                   | 60                | 6.5-7.5                     | 6                                  | SF08060SV                                     | SF08060MV |
| 8                   | 80                | 6.5-7.5                     | 6                                  | SF08080SV                                     | SF08080MV |
| 8                   | 100               | 6.5-7.5                     | 6                                  | SF08100SV                                     | SF08100MV |
| 8                   | 120               | 6.5-7.5                     | 6                                  | SF08120SV                                     | SF08120MV |
| 8                   | 150               | 6.5-7.5                     | 6                                  |                                               | SF08150MV |
| 8                   | 200               | 6.5-7.5                     | 6                                  | SF08200SV                                     | SF08200MV |



For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and S.M.A.R.T. are trademarks of Cordis and may be registered in the US and/or in other countries.  
© 2023 Cordis. All Rights Reserved. 100551514-4 01/2023