

# Cordis **STABILIZER**™ Support Guidewire and Cordis **STABILIZER**™ Extra Support Guidewire 0.014" Peripheral Steerable Guidewires



Providing All the Support  
You Need for 0.014"  
Peripheral Interventions



# Cordis **STABILIZER™** Support Guidewire and Cordis **STABILIZER™** Extra Support Guidewire

## 0.014" Peripheral Steerable Guidewires

### Versatility

Broad transitions and excellent support for device deliverability and compatibility with OUTBACK™ LTD™ Re-Entry Catheter

### Enhanced Lubricity

Proximal DURAGLIDE™ PTFE Spray and distal PTFE Sleeve for exceptional surface lubricity

### Atraumatic Tip

SuperSoft tip flexibility provides atraumatic vessel entry and navigation

## Ordering Information

### Cordis **STABILIZER™** Support Guidewire

Catalog Numbers	Diameter (inches)	Tip Flexibility	Length (cm)	Tip Shape
507 - 180S	0.014	SuperSoft	180	Straight
507 - 300S	0.014	SuperSoft	300	Straight

### Cordis **STABILIZER™** Extra Support Guidewire

Catalog Numbers	Diameter (inches)	Tip Flexibility	Length (cm)	Tip Shape
527 - 180E	0.014	SuperSoft	180	Straight
527 - 300E	0.014	SuperSoft	300	Straight

### Compatible with **OUTBACK™ LTD™** Re-Entry Catheter

Codes	Length (cm)	Sheath (F)	-	-
OTB 42120	120	6F Sheath	-	-

### Cordis **STABILIZER™** Support Guidewire

Support for the delivery of advanced peripheral devices in complex and tortuous anatomy



### Cordis **STABILIZER™** Extra-Support Guidewire

When a higher level of support is needed for the most challenging peripheral interventions

