


Cordis EMERALD™

Diagnostic Guidewires



PTFE pre-coated
guidewire

Rosen Heavy Duty and Amplatz
Super Stiff EMERALD™ Guidewires

Rosen Heavy Duty and Amplatz Super Stiff PTFE-Coated EMERALD™ Guidewires*

EMERALD™ Amplatz Super Stiff Guidewires

Diameter (inches)	Length (cm)	J-Radius (mm)	Catalogue Number
PTFE-Coated STRAIGHT TIP			
.035 Amplatz	150	–	502-726
.035 Amplatz	180	–	502-728

PTFE-Coated J-CURVE

.035 Amplatz	150	3	502-731
.035 Amplatz	180	3	502-733
.035 Amplatz	260	3	502-735


- For percutaneous entry and guidance of angiographic catheters
- Pre-coated PTFE
- 5 units per package

EMERALD™ Rosen Heavy Duty Guidewires

Diameter (inches)	Length (cm)	J-Radius (mm)	Catalogue Number
PTFE-Coated J-CURVE			
.035 Rosen	180	1.5	502-717

- Atraumatic end
- For percutaneous entry and high support guidance of angiographic catheters
- 5 units per package

* More configurations available. Please consult your Cordis Catalogue.



For more information please visit www.cardinalhealth.co.uk



For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO and EMERALD are trademarks of Cardinal Health and may be registered in the US and/or in other countries.

© 2021 Cardinal Health. All Rights Reserved. 100552138-2 02/2021

