


# NYLEX<sup>®</sup>

4F & 5F



## High Quality non-braided Flush Catheters



- Kink resistance
- Ease of introduction
- High flow
- Excellent flexibility
- Superb shape retention
- Guidewire acceptance .035"

# NYLEX™ available in 4F & 5F

## Performance Makes the Difference

Shape	Side Holes	Length (cm)	Flow Rates		Catalogue Number	
			4F	5F	4F	5F
Pigtail	8	65	17.41 ml/s	33 ml/s	526-410	526-510
Pigtail	8	90	15.87 ml/s	28 ml/s	526-411	526-511
Pigtail	8	100	15.80 ml/s	26 ml/s	526-412	526-512
Pigtail	8	110	15.10 ml/s	24 ml/s	526-413	526-513
Pigtail	8	65	17.41 ml/s	33 ml/s	526-430	526-530
Pigtail	8	90	15.87 ml/s	28 ml/s	526-431	–
Straight	8	65	17.41 ml/s	33 ml/s	526-420	526-520
Straight	8	90	15.87 ml/s	28 ml/s	526-421	526-521
Straight	8	100	15.80 ml/s	26 ml/s	526-422	526-522
Straight	8	110	15.10 ml/s	24 ml/s	526-423	526-523
Universal Flush	8	65	17.41 ml/s	33 ml/s	526-440	526-540
Universal Flush	8	100	15.80 ml/s	26 ml/s	526-442	–
D.S.A. set*, Pigtail	8	65	17.41 ml/s	33 ml/s	503-410	503-510
D.S.A. set*, straight	8	65	17.41 ml/s	33 ml/s	503-420	503-520

\* including .035" guidewire and 19 G needle

*Inner lumen • 4F from hub to tip .040" • 5F .050" body & .038" tip  
• Max. injection pressure 1050 psi*

### Outstanding Tip Characteristics

- radiopacity to facilitate positioning, save time and contrast medium
- tapered tip (5F) to facilitate skin penetration

### Kink Resistance, Excellent Flexibility

- pigtail easily straightens for smooth placement
- ease of introduction (no need for a vessel dilator)

### Superb Shape Retention

- pigtail accurately regains its original shape to optimize dispersion of contrast dye and stability during injection

**Important Information:** Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings, and precautions. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. For healthcare professionals only.

© Cardinal Health 2016. All right reserved. CORDIS, the Cordis LOGO, NYLEX™ are trademarks or registered trademarks of Cardinal Health.

• EU507-1 10/16

**Cordis**  
A Cardinal Health company