

PALMAZ GENESIS™ Peripheral Stent on OPTA™ PRO .035" Delivery System (Large)

.035" Over-the-wire System PALMAZ GENESIS™ on OPTA™ PRO (Large)

- Recommended inflation pressure : 8 atm



designed for life

Ordering Information

Balloon diameter x length (mm x mm)	Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Sheath Compatibility (F)	Guiding Compatibility (F)	Catalog Number Length of Delivery System	
					80 cm	135 cm
5 x 40	39	39	6	8 (.088")	PG3950PPS	PG3950PPX
5 x 60	59	59	6	8 (.088")	PG5950PPS	PG5950PPX
6 x 30	29	27	6	8 (.088")	PG2960PPS	PG2960PPX
6 x 40	39	38	6	8 (.088")	PG3960PPS	PG3960PPX
6 x 60	59	59	6	8 (.088")	PG5960PPS	PG5960PPX
6 x 80	79	79	6	8 (.088")	PG7960PPS	PG7960PPX
7 x 30	29	26	6	8 (.088")	PG2970PPS	PG2970PPX
7 x 60	59	57	6.5	8 (.088")	PG5970PPS	PG5970PPX
7 x 80	79	78	6.5	8 (.088")	PG7970PPS	PG7970PPX
8 x 40	39	36	7	9 (.098")	PG3980PPS	PG3980PPX
8 x 60	59	56	7	9 (.098")	PG5980PPS	PG5980PPX
8 x 80	79	76	7	9 (.098")	PG7980PPS	PG7980PPX
9 x 20	19	17	7	9 (.098")	PG1990PPS	PG1990PPX
9 x 26	25	22	7	9 (.098")	PG2590PPS	PG2590PPX
9 x 30	29	25	7	9 (.098")	PG2990PPS	PG2990PPX
9 x 40	39	35	7	9 (.098")	PG3990PPS	PG3990PPX
9 x 60	59	53	7	9 (.098")	PG5990PPS	PG5990PPX
9 x 80	79	75	7	9 (.098")	PG7990PPS	PG7990PPX
10 x 20	19	16	7	9 (.098")	PG1910PPS	PG1910PPX
10 x 26	25	22	7	9 (.098")		PG2510PPX
10 x 30	29	24	7	9 (.098")	PG2910PPS	PG2910PPX
10 x 40	39	32	7	9 (.098")	PG3910PPS	PG3910PPX
10 x 60	59	53	7	9 (.098")	PG5910PPS	PG5910PPX
10 x 80	79	73	7	9 (.098")	PG7910PPS	

Complementary Cordis Products


The PALMAZ GENESIS™ Stent on OPTA™ PRO .035" Delivery System can be used with other Cordis product family. The products listed hereunder are only a selection. Please contact your local sales representative for further information on product compatibilities.

The PALMAZ GENESIS™ on OPTA™ Pro Stent can be used with other Cordis product family. The products listed hereunder are only a selection.

Ordering Information

.035" STORQ™

Steerable Guidewires



Diameter (in)	Tip Flexibility	Tip Shape	Length (cm)	Catalogue Number
.035	Standard	Straight	180	503-456
.035	Standard	Angled	180	503-456J
.035	Standard	Modified J	180	503-456MJ
.035	Standard	Straight	300	503-456X
.035	Standard	Angled	300	503-456Y
.035	Standard	Modified J	300	503-456MY

BRITE TIP™

Sheath Introducers



Cannula usable length (cm)	Catalogue Number 6F	Catalogue Number 7F	Catalogue Number 8F
11	401-611 M	401-711 M	401-811 M
23	401-623 M	401-723 M	401-823 M
35	401-635 M	401-735 M	401-835 M
45	401-645 M	401-745 M	401-845 M
55	401-655 M	401-755 M	401-855 M
90	401-690 M	401-790 M	

MRI Compatibility (Please refer to the IFU-section IX)



For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO, BRITE TIP, STORQ, OPTA PRO and PALMAS GENESIS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100549930-2 01/2021

For more information please visit www.cardinalhealth.co.uk



Cordis
A Cardinal Health company