

PALMAZ GENESIS™ on OPTA™ Pro Stent

(Medium)

.035" Over-the-wire System

PALMAZ GENESIS™ on OPTA™ Pro (Medium)

- Recommended inflation pressure: 8 atm



designed for life

Ordering Information

Balloon diameter x length (mm x mm)	Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Sheath Compatibility (F)	Guiding Compatibility 8F / .088"	Catalog Number	
					Length of Delivery System 80 cm	135 cm
4 x 20	18	17	6	8	PG1840PPS	
4 x 25	24	23	6	8	PG2440PPS	
5 x 17	15	14	6	8	PG1550PPS	
5 x 20	18	17	6	8	PG1850PPS	PG1850PPX
5 x 25	24	23	6	8	PG2450PPS	
6 x 17	15	14	6	8	PG1560PPS	
6 x 20	18	17	6	8	PG1860PPS	PG1860PPX
6 x 25	24	23	6	8	PG2460PPS	PG2460PPX
7 x 15	12	11	6	8	PG1270PPS	PG1270PPX
7 x 17	15	14	6	8	PG1570PPS	
7 x 20	18	17	6	8	PG1870PPS	
7 x 25	24	23	6	8	PG2470PPS	PG2470PPX
8 x 15	12	11	6	8	PG1280PPS	PG1280PPX
8 x 17	15	14	6	8	PG1580PPS	PG1580PPX
8 x 20	18	17	6	8	PG1880PPS	PG1880PPX
8 x 25	24	23	6	8	PG2480PPS	

Complementary Cordis Products


The PALMAZ GENESIS™ on OPTA™ Pro Stent can be used with other Cordis product family. The products listed hereunder are only a selection. please contact your local sales representative for further information on product compatibilities.

The PALMAZ GENESIS™ on OPTA™ Pro Stent can be used with other Cordis product family. The products listed hereunder are only a selection.

Ordering Information

JINDO™


Steerable Guidewires



Diameter (in)	Taper	Length (cm)	Catalogue Numbers
.035	Short	180	503-452
.035	Short	300	503-453
.035	Intermediate	300	503-553

VISTA BRITE TIP™ IG

Introducer Guides



Catheter Shape	Length (cm)	OD (F)	ID	Catalogue Number
Renal curve (small)	55	8	.088"	403-8553R
Renal curve (long)	55	8	.088"	403-8553A

