



MORE CASES
COVERED

2-8MM DIAMETERS

20-300MM LENGTHS

MORE .018" RX BALLOON SIZES

SIZES THAT TREAT LESIONS AND LOCATIONS



LARGE LESIONS

- DIAMETERS 6-8mm
- LENGTHS 20-100mm
- 4-6 FRENCH FIT

HIGHLY DIFFUSE LESIONS

- DIAMETERS 4-6mm
- LENGTHS UP TO **300mm**
- 4-6 FRENCH FIT

SMALL LESIONS


- DIAMETERS 2-4mm
- LENGTHS UP TO **300mm**
- SHAFT LENGTHS UP TO 155cm
- 4 FRENCH FIT


DESIGNED WITH EXCELLENT CROSSABILITY

Molded tip design and low-profile body

To enhance crossability in
challenging lesions

Dual layer hydrophilic coated balloon and outer shaft





Cordis SABERX™ Catheter offers:

- Designed proximal shaft with a tapered core wire
- 40cm OTW length allows the RX port remain in the long sheath introducer/guiding catheter in many lower limb cases

DESIGNED FOR **EXCELLENT CROSSABILITY**

- Dual-layer hydrophilic coating
- Molded soft tip with outstanding transition to wire
- Low profile design

ENHANCED **PERFORMANCE FEATURES**

- Designed guidewire exit transition area for higher kink resistance
- RBP up to 18 atm
- DURALYN® Balloon Material

2-8MM DIAMETERS

20-300MM LENGTHS

CORDIS SABERX™ PTA DILATATION CATHETER

ORDERING INFORMATION

Product Description

Type	Low profile, high pressure, rapid-exchange PTA balloon, intended to dilate stenoses in iliac, femoral, ilio-femoral, popliteal, infra popliteal and renal arteries.
Shaft Diameter	3.9F and 3.3F
Guidewire	0.018"
Usable Lengths	155cm
Coating	Dual layer, hydrophilic coating
Features	MORE 0.018" Rapid Exchange sizes.

Ordering Information

Balloon Diameter, mm	Balloon Length, mm	Rated Burst Pressure, atm	Sheath Fit, (F)	Shaft Length (155cm)
2.0	20	18	4	51002002L
2.0	30	18	4	51002003L
2.0	40	18	4	51002004L
2.0	60	18	4	51002006L
2.0	80	18	4	51002008L
2.0	100	18	4	51002010L
2.0	150	18	4	51002015L
2.0	200	18	4	51002020L
2.0	250	18	4	51002025L
2.0	300	18	4	51002030L
2.5	20	18	4	51002502L
2.5	30	18	4	51002503L
2.5	40	18	4	51002504L
2.5	60	18	4	51002506L
2.5	80	18	4	51002508L
2.5	100	18	4	51002510L
2.5	150	18	4	51002515L
2.5	200	18	4	51002520L
2.5	250	18	4	51002525L
2.5	300	18	4	51002530L
3.0	20	18	4	51003002L
3.0	30	18	4	51003003L
3.0	40	18	4	51003004L
3.0	60	18	4	51003006L
3.0	80	18	4	51003008L
3.0	100	18	4	51003010L
3.0	150	14	4	51003015L
3.0	200	14	4	51003020L
3.0	250	14	4	51003025L
3.0	300	14	4	51003030L
3.5	20	18	4	51003502L
3.5	30	18	4	51003503L
3.5	40	18	4	51003504L
3.5	60	18	4	51003506L
3.5	80	18	4	51003508L
3.5	100	18	4	51003510L
3.5	150	14	4	51003515L
3.5	200	14	4	51003520L
3.5	250	14	4	51003525L
3.5	300	14	4	51003530L
4.0	20	18	4	51004002L

Balloon

Material	DURALYN®
Markers	Proximal and distal marker bands (dual proximal marker bands for balloon length ≥ 100 mm)
Nominal	6 to 8 atm
RBP	12 to 18 atm
Shoulder Angle	Smooth balloon shoulders
Balloon Lengths	20mm, 30mm, 40mm, 60mm, 80mm, 100mm, 150mm, 200mm, 250mm, 300mm
Balloon Diameters	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 5mm, 6mm, 7mm, 8mm

Balloon Diameter, mm	Balloon Length, mm	Rated Burst Pressure, atm	Sheath Fit, (F)	Shaft Length (155cm)
4.0	30	18	4	51004003L
4.0	40	18	4	51004004L
4.0	60	18	4	51004006L
4.0	80	18	4	51004008L
4.0	100	18	4	51004010L
4.0	150	14	4	51004015L
4.0	200	14	4	51004020L
4.0	250	14	4	51004025L
4.0	300	14	4	51004030L
5.0	20	16	4	51005002L
5.0	30	16	4	51005003L
5.0	40	16	4	51005004L
5.0	60	16	4	51005006L
5.0	80	16	4	51005008L
5.0	100	16	4	51005010L
5.0	150	14	5	51005015L
5.0	200	14	5	51005020L
5.0	250	14	5	51005025L
5.0	300	14	5	51005030L
6.0	20	16	4	51006002L
6.0	30	16	4	51006003L
6.0	40	16	4	51006004L
6.0	60	16	4	51006006L
6.0	80	16	5	51006008L
6.0	100	16	5	51006010L
6.0	150	14	5	51006015L
6.0	200	14	5	51006020L
6.0	250	14	6	51006025L
6.0	300	14	6	51006030L
7.0	20	14	4	51007002L
7.0	30	14	4	51007003L
7.0	40	14	4	51007004L
7.0	60	14	5	51007006L
7.0	80	14	5	51007008L
7.0	100	14	5	51007010L
8.0	20	12	4	51008002L
8.0	30	12	4	51008003L
8.0	40	12	4	51008004L
8.0	60	12	5	51008006L
8.0	80	12	5	51008008L
8.0	100	12	5	51008010L

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and SABERX are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. © 2020 Cardinal Health. All Rights Reserved. 100557683 03/2020

