

# SUPER TORQUE<sup>®</sup> MB

5F Marker Band Flush Catheters



The Gauge  
for Success




# SUPER TORQUE® MB

## 5F Marker Band Flush Catheters

Matching the exacting standards of Cordis quality, the SUPER TORQUE® MB markerband catheter adds the possibility of performing precise length measurements to the angiographic features of the SUPER TORQUE® MB diagnostic catheters.

The excellent visibility of the gold alloy marker bands provide easy location information of the catheter in the vasculature.

The braided body provides you with a highly manoeuvrable device that is demanded by this type of challenging procedure.



### Ordering Information

Lengths (cm)	Size (F)	2 gold-alloy marker bands		10 gold-alloy marker bands		20 gold-alloy marker bands	
		SH*	Cat. N°	SH*	Cat. N°	SH*	Cat. N°
65	5	–		–		8	532-598C
70	5	6	532-598D**	–		6	
110	5	6		6	532-598A	6	532-598B

\*Side Holes

\*\* 2 cm spacing per beginning of each marker