

PALMAZ GENESIS™ Peripheral Stent on SLALOM™ .018" Delivery System



.018" Over-the-wire System PALMAZ GENESIS™ on SLALOM™

Ordering Information

Balloon Diameter x Length (mm x mm)	Unexpanded Stent Length (mm)	Expanded Stent length (mm)	Rated Burst Pressure (atm)	Nominal Pressure (atm)	Crossing Profile (inch.)	Sheath fit (F)	Guiding Catheter fit (F)	Catalog Number	
								Length of Delivery System 80 cm	Length of Delivery System 135 cm
3 x 20	18	18	14	12	.055"	5	6 (.067")		PG1830PSX
3 X 25	24	24	14	12	.055"	5	6 (.067")		PG2430PSX
7 x 40	39	37	12	10	.074"	6	7 (.078")		PG3970PSX
8 x 20	18	16	10	8	.080"	6	8 (.088")		PG1880PSX
8 x 25	24	23	10	8	.080"	6	8 (.088")		PG2480PSX
8 x 40	39	34	10	8	.080"	6	8 (.088")	PG3980PSS	PG3980PSX

Complementary Cordis Products

The PALMAZ GENESIS™ on SLALOM™ Stent can be used with other CORDIS™ products family. The products listed hereunder are only a selection. Please contact your local sales representative for further information on product compatibilities.

Compatible guidewires diameter are listed on the label.


Ordering Information

.018" SV Wire Steerable Guidewires



Guidewire	Taper length (cm)	Length (cm)	Tip	Catalogue Number
SV-5	5	180	Straight	503-558
SV-5	5	300	Straight	503-558X

VISTA BRITE TIP™ Guiding Catheters



Catheter Shape	Length (cm)	Catalogue Number 6F" (.067")
Renal Curve (small)	55	670-210-55
Renal Curve (long)	55	670-212-55

