



Cordis **ELITECROSS**TM

Support Catheter



**Power Steering for
More Powerful
Crossing.**

Providing a unique combination of support and shape options to enable a tailored approach to tough lesions.


- UNPARALLELED SUPPORT*
- UNPARALLELED PUSHABILITY*
- ENHANCED STEERABILITY*
- EXTRA LARGE INNER LUMEN can accommodate 0.014" support catheter systems
- Tapered tip and lubricious hydrophilic coating enhance deliverability

Portfolio:


- 0.035" guidewire compatibility
- 82 cm and 132 cm lengths
- XS Straight and XS Angled versions
- DX and IX device compatible

Configuration Comparison:

ELITECROSS™ Support Catheter - Straight EXTRA SUPPORT



ELITECROSS™ Support Catheter - Angled EXTRA SUPPORT



ELITECROSS™ SUPPORT CATHETER

Catalog Code	Description	Shape	Length
MGC39090X	ELITECROSS™ SUPPORT CATHETER - EXTRA SUPPORT	STRAIGHT	82 cm
MGC39090A	ELITECROSS™ SUPPORT CATHETER - EXTRA SUPPORT	ANGLED	82 cm
MGC39140X	ELITECROSS™ SUPPORT CATHETER - EXTRA SUPPORT	STRAIGHT	132 cm
MGC39140A	ELITECROSS™ SUPPORT CATHETER - EXTRA SUPPORT	ANGLED	132 cm

* Cordis internal testing reports

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

For healthcare professionals only. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. Contact your Cordis sales representative for availability and ordering. Important information: Prior to use, refer to the Instructions for Use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings, and precautions. © 2016 Cardinal Health. All rights reserved. CORDIS, the Cordis LOGO, and ELITECROSS are trademarks or registered trademarks of Cordis Corporation. 155-8839 23574 1/16