



POWERFLEX[®] PRO

0.035" PTA Dilatation Catheter

Advanced crossability and
versatility to treat routine, or
challenging cases.



POWERFLEX® PRO

0.035" PTA Dilatation Catheter

ACCESS


- Low Profile:
 - » Wide array of 5F compatible sizes
 - » Excellent trackability and pushability to facilitate lesion access.
- Kink resistant shaft:
 - » To facilitate pushability when encountering tortuosity.
 - » To help avoid procedure delays.

ADVANCED CROSSABILITY

- Low profile and atraumatic flexible tip to facilitate lesion access.
- MDX Silicone coated Distal (30cm) shaft and guidewire lumen to facilitate smooth delivery and wire handling.
- DURALYN™ combines flexibility and strength for improved crossability and dilation of long lesions.

VERSATILITY

- Powerful balloon inflation with Rated Burst Pressures up to 18 ATM.
- DURALYN™ High Performance Balloon Material provides controlled compliance for predictable balloon sizing.
- Rapid inflation and deflation time to maximize procedural efficiency.
- Radiopaque markers facilitate placement of the POWERFLEX® PRO PTA Catheter for effective dilatation.



FEATURES

- Dual lumen shaft
- Atraumatic tip
- RBP up to 18 ATM
- Extended product range
- One piece hub
- Inner and outer (30cm distal) MDX coating

ORDERING INFORMATION

Packaging: 1 unit/Box

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (80cm)	Shaft Length (135cm)
3	20	18	5	4400302S	4400302X
3	40	18	5	4400304S	4400304X
3	60	18	5	4400306S	4400306X
3	80	18	5	4400308S	4400308X
3	100	18	5	4400310S	4400310X
3	120	18	5	4400312S	4400312X
3	150	18	5	4400315S	4400315X
3	220	18	5	4400322S	4400322X
4	20	18	5	4400402S	4400402X
4	40	18	5	4400404S	4400404X
4	60	18	5	4400406S	4400406X
4	80	18	5	4400408S	4400408X
4	100	18	5	4400410S	4400410X
4	120	18	5	4400412S	4400412X

POWERFLEX® PRO

0.035" PTA Dilatation Catheter

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (80cm)	Shaft Length (135cm)
4	150	18	5	4400415S	4400415X
4	220	18	5	4400422S	4400422X
5	20	15	5	4400502S	4400502X
5	40	15	5	4400504S	4400504X
5	60	15	5	4400506S	4400506X
5	80	15	5	4400508S	4400508X
5	100	15	5	4400510S	4400510X
5	120	15	5	4400512S	4400512X
5	150	15	5	4400515S	4400515X
5	220	15	5	4400522S	4400522X
6	20	15	5	4400602S	4400602X
6	40	15	5	4400604S	4400604X
6	60	15	5	4400606S	4400606X
6	80	15	5	4400608S	4400608X
6	100	15	5	4400610S	4400610X
6	120	15	5	4400612S	4400612X
6	150	15	5	4400615S	4400615X
6	220	15	5	4400622S	4400622X
7	20	15	5	4400702S	4400702X
7	40	15	5	4400704S	4400704X
7	60	12	5	4400706S	4400706X
7	80	12	5	4400708S	4400708X
7	100	12	6	4400710S	4400710X
8	20	12	5	4400802S	4400802X
8	40	12	5	4400804S	4400804X
8	60	12	5	4400806S	4400806X
8	80	12	5	4400808S	4400808X
8	100	12	5	4400810S	4400810X
9	20	12	6	4400902S	4400902X
9	40	12	6	4400904S	4400904X
9	60	12	6	4400906S	4400906X
9	80	12	7	4400908S	4400908X
9	100	12	7	4400910S	4400910X
10	20	12	6	4401002S	4401002X
10	40	12	6	4401004S	4401004X
10	60	12	7	4401006S	4401006X
10	80	12	7	4401008S	4401008X
10	100	12	7	4401010S	4401010X
12	20	8	7	4401202S	4401202X
12	40	8	7	4401204S	4401204X
12	60	8	7	4401206S	4401206X

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions.

CORDIS, Cordis LOGO, POWERFLEX, and DURALYN are registered trademarks of Cardinal Health and may be registered in the US and/or in other countries.

© 2021 Cardinal Health. All Rights Reserved. 100586670 07/2021

Cordis
A Cardinal Health company