


RADIAL-SPECIFIC DESIGN. PATIENT-FOCUSED BENEFITS.

SABERX RADIANTZ™ Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter

Clinically comprehensive PTA-size matrix in an .018" rapid-exchange catheter that can be delivered radially¹




Designed for Radial¹

190 cm shaft length facilitates treatment of iliac and SFA lesions via a radial approach

Low-profile 4Fr sheath compatibility for most sizes reduces the potential for access site complications

Rapid Exchange delivery system in an .018" platform enables convenient device exchange and optimized wire control

Complete iliac matrix including 9-10 mm diameter balloons facilitates pre/post-dilatation complementing S.M.A.R.T. RADIANTZ™ Vascular Stent System



Engineered for Exceptional Deliverability

Low-crossing profile and atraumatic tapered tip are designed to enhance crossability

Dual-layer hydrophilic coated balloon and outer shaft reduces friction for enhanced trackability and crossability

Breakthrough Tack Seal shaft design enhances pushability and reduces risk of kinking



Streamlined for Procedural Efficiency

300 mm balloon lengths enable treatment of long lesions with fewer inflations and may potentially reduce procedural time

Trusted DURALYN™ Material provides controlled compliance and smooth shoulders resulting in greater strength


High RBP (up to 18 atm) enables effective dilatation of challenging lesions¹

¹ Cordis 2020 Data on File

SABERX RADIANTZ™ Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter

SABERX RADIANTZ™ Catheter provides superior crossability over competing .018" radial PTA catheters (p < .05)

Extensive size matrix facilitates broad treatment of lesions and may reduce procedural time




*4 x 100 mm balloon tested in a Radial Peripheral Crossing Model

Balloon Lengths (mm)

	20	30	40	60	80	100	150	200	250	300
2.0	18	18	18	18	18	18	18	18	18	18
2.5	18	18	18	18	18	18	18	18	18	18
3.0	18	18	18	18	18	18	14	14	14	14
3.5	18	18	18	18	18	18	14	14	14	14
4.0	18	18	18	18	18	18	14	14	14	14
5.0	16	16	16	16	16	16	14	14	14	14
6.0	16	16	16	16	16	16	14	14	14	14
7.0	14	14	14	14	14	14				
8.0	12	12	12	12	12	12				
9.0	10	10	10	10	10	10				
10.0	10	10	10	10	10	10				

Legend: 4F (purple), 5F (grey), 6F (green)
Values in table represent RBP

Product Overview



RADIAL-SPECIFIC DESIGN. PATIENT-FOCUSED BENEFITS.

Learn more at www.cordis.com or contact your local Cordis sales rep or customer service at 800.327.7714.